

Bio-hazard Kits - Infection Control

“ Effective handling of infectious bodily fluids, with Steroplast’s Bio-hazard Kits ”



Contents - Ref BIO2

| | |
|---------------------|----------------------------|
| Gloves & apron: | protects against infection |
| Absorbent granules: | solidifies the spill |
| Scoop & scraper: | aides removal of the spill |

| | |
|----------------|-----------------------------|
| Waste bag: | for all debris |
| Steroclenz 50: | for surface decontamination |
| Wipe: | to clean hands |



Tub of super absorbent powder 100g



Open the camera app on your smartphone device and hover over the QR code to find out more

The sachet/tub of super absorbent powder is poured onto the spill **which immediately starts to turn the liquid, into a solid for easy removal.**

The yellow scoop and scraper are provided in the kit to make cleaning up the solidifying spill quick and easy. The waste can now be placed in the yellow clinical waste bag provided.

Now disinfect the area with the disinfecting spray provided. Place the used scoop, scraper and gloves in the yellow Bio-hazard bag. Dispose of yellow Bio-hazard bag according to your local authority/council guidelines.

Faced with highly infectious, and potentially life threatening bodily fluids to clean up, you'll need a Steroplast Bio-hazard kit for your protection.

Bodily fluids including blood, vomit, urine, mucus, faeces, semen either in isolation, or on objects such as needles, syringes, broken glass and knives can all be safely disposed of using the following kits.

How to use the Steroplast Bio-hazard kit to clean up bodily fluids, quickly and safely:

- 1** Using the sachet of super absorbent powder/granules, pour this onto the spillage area. This will immediately take effect and the liquid will begin to solidify. This makes it much easier to remove.
- 2** The yellow scoop and scraper provided in the kit make cleaning up the solidifying spill quick and easy. The gathered waste can now be placed into the yellow clinical waste bag provided in the kit.
- 3** Now disinfect the area with the spray provided. Place the used scoop, scraper and gloves into the yellow biohazard bag provided.



Ref BIO2

| Description | Item | Ref |
|--------------------------------------|------|--------|
| Clean up kit - 5 applications | 1 | BIO5 |
| Clean up kit - 3 applications | 1 | BIO3 |
| Clean up pack - 1 application | 1 | BIO2 |
| Clean up pack refill | 1 | BIO1-R |
| Absorbent powder - 100g tub | 1 | 8213 |

For safe and effective handling, containment, disinfection and disposal of infectious microorganisms contained in blood, vomit, urine, mucus, faeces and semen.

The kits meets all Health & Safety guidelines.



Ref BIO3

Designed for up to 3 separate incidents