

iPAD SP1 Defibrillator



The iPAD SP1 can save a life, with its innovative features and advanced technology

Every second counts

The iPAD uses the latest in defibrillation technology to ensure that the device works as quickly as possible. After a cycle of 2 minutes of CPR, the device can perform analysis and be ready for shock delivery in less than 5 seconds.

Simple to operate and maintain

There is guidance at every step of the way with clear illuminated illustrations, calm voice prompts and large responsive buttons. The iPAD SP1 intelligently guides a rescuer at the pace at which they work.

Simple switch selection from adult to child mode

Just using a simple switch system the same set of electrode pads can be

used for both adults and children.

There are no additional pads required or setting changes, saving precious time and expense.

Background noise is detected and controlled

The iPAD SP1 has an effective, ambient noise detector with voice prompts adjustable up to 90 decibels. This ensures that the machine's important verbal guidance will always be heard above high levels of background noise.

Integrates seamlessly with the Ambulance Service

The iPAD SP1 is totally compatible with front line emergency vehicles. It contains plug-in pad adapters to ensure a seamless handover of a patient to advanced care.

Extremely well made – solid and reliable

The iPAD SP1 has been manufactured to a very high standard. It has a robust one piece design and there is nothing to break or lose. It can be relied on to perform and endure emergency usage time after time.

To order call:

0161 902 3030

for more info email:

marketing@steroplast.co.uk

or visit:

www.steroplast.co.uk

 To order call **0161 902 3030**

 email **marketing@steroplast.co.uk**

 or visit **www.steroplast.co.uk**



Setting Standards Across Seven Sectors

Emergency

Hospital

Sport

Veterinary

Medical

First Aid

Retail