

Diagnostics

Here is a reflection of the diagnostic equipment available at Steroplast, visit our website or see QR code for full list.



Open the camera app on your smartphone device and hover over the QR code to find out more

Forehead thermometer

Available in 25's and 50's

The unbreakable and easy to use forehead thermometer gives fast and accurate readings.

Advantages, features & benefits

- Completely safe
- The pocket-sized product can be taken anywhere
- Fast and simple to use
- Temperature range: 35 - 40 Celsius



Ref 1811/25
Ref 1811/50

- Thermometer dimensions: 89mm x 19mm
- Latex-Free
- Hospital and home use
- The fast and simple temperature monitoring system can be used as pre and post-operation surveillance.

Digital ear thermometer

Compact, ergonomic digital thermometer with a 1 second measurement response time.

- Holds a capacity of 10 memory records
- Auto shut off
- Approx. 4,000 uses between batteries
- Fast and accurate



Ref 1128



Ref 1802

Braun ThermoScan ear thermometers

Ear thermometer measures the infrared heat generated by the eardrum and surrounding tissue.

- Dual scale (celsius or fahrenheit)
- Ultra soft and gentle tip is warmed for extra comfort



Ref 1817



Ref 1799

Digital thermometer

A safe, accurate and simple way to check and monitor your body temperature.

- Accurate and easy to read display
- Simple to operate
- No glass or mercury
- Provides fast and accurate body temperature reading
- Large screen for fast and easy reading of results
- Automatic shut-off function preserves battery life



Ref 1818

Berrcom non-contact infrared thermometer

This forehead thermometer is designed to take a temperature with a 3-5cm distance from the forehead. It uses infrared technology to measure temperature and is quick and easy to use.

- Non-contact for hygienic, non-intrusive use
- Audible fever alarm will alert you to a temperature of over 38.6°C
- Easy to read, backlit LCD screen
- Temperature displayed in either Fahrenheit or Celsius
- Body measurement range: 25°C - 42.9°C
- Auto power off after 30 seconds



Ref 1796

Twin tube Sprague Rappaport stethoscope

Traditional design stethoscope with two tubes connecting the chest piece to the head frame.

- Two tubes connecting the chest piece to the head frame
- Supplied with interchangeable chest pieces for adult and paediatric use
- Available in 6 additional colours Red/Purple/Pink/Orange/Blue/Green/Black

Ref WSTE0116/RED



Cardiology stethoscope

- Stainless steel double head chest piece with a deep cone-shaped bell gives excellent acoustic sensitivity
- The thick tube offers resistance to ambient interference
- Binaural with reinforced stainless steel inner-spring for durability and portability
- 3-year warranty



Ref 1063

Deluxe cardiology stethoscope

- Excellent sound conduction
- Zinc alloy chest piece
- Flexible binaural for optimal user comfort
- Durable soft silicone ear tips
- Supplied with spare ear tips
- 3-year warranty



Ref 1118

More Stethoscopes available on our website www.steroplast.co.uk

Doctors surgery and hospitals within a 20 mile radius from our office, can have products delivered within 2 hours

Creative medical fingertip pulse oximeter

The tough, reliable and accurate PC-60B1 EMS Professional Model.

- 3 Year drop warranty - Super tough EMS/Military model
- Automatic ON/OFF for ease of use
- 4 Direction dual colour LCD display
- Plethysmographic waveform allows optimum finger placement
- Large backlit display of SpO2, PR, PI and pulse bar
- Supplied with carry pouch, lanyard and batteries



Ref WSTE0077

Omron M6 AC LED blood pressure monitor

Clinically validated monitor incorporates the most advanced Omron sensors and indicators for accurate and reliable blood pressure readings.

- High performance, quality digital blood pressure monitor
- 100 reading, 2 user memory
- Body movement detector
- Hypertension tracker
- Clinical and Diabetic validation
- Stroke prevention



Monitor Cuff - Large (Ref WSTE0111)

Freestyle Optium Neo blood glucose meter

Designed to make life simpler for people who use insulin.

- Accurate blood glucose and blood ketone testing
- Insulin logging helps track insulin doses
- Fast, five second test time
- Notification when results need attention
- Clear and sharp screen
- Matte screen reads like paper, no glare
- Icon driven display

Ref FFR2029

