

Steroply

Conforming Bandage

A high stretch, lightweight bandage for holding dressings in place.

- Used in NHS hospitals
- Soft edge for comfort
- Highly conformable, fits snugly anywhere on the body
- A nylon/viscose construction creates fray resistance
- Light and soft, permitting air circulation
- Individually wrapped
- Choice of sizes

Steroply is renowned for its quality



Description	Size	Pack	Ref
Steroply	5cm x 4m	1 x 12	3051
Steroply	7.5cm x 4m	1 x 12	3052
Steroply	10cm x 4m	1 x 12	3053
Steroply	15cm x 4m	1 x 12	3054

NHS Approved


Sterocrepe

Crepe Bandage

The NHS approved crepe bandage, ideal for giving light support to sprains and strains.

- Unrivalled quality
- Superb elasticity through "Z" and "S" twist yarns
- Can also be used for dressing retention
- Soft and supple
- Washable and therefore re-usable
- Individually wrapped

Used extensively in Health & Beauty for body wraps





Description	Size	Pack	Ref
Sterocrepe	5cm x 4.5m	1 x 12	3056
Sterocrepe	7.5cm x 4.5m	1 x 12	3057
Sterocrepe	10cm x 4.5m	1 x 12	3058
Sterocrepe	15cm x 4.5m	1 x 12	3059

Wide-open-weave (WOW) Bandage

The best selling bandages with a wide range of applications

- Loose weave construction allows wounds to breathe
- 100% cotton - Non elastic - available in two widths



Description	Size	Pack	Ref
Wow - bandage	5cm x 5m	1 x 12	3048
Wow - bandage	7.5cm x 5m	1 x 12	3049

Triangular Bandages

The multifunctional bandage can support an arm in a sling or as additional support to help hold a splint in place

- Available in a non-woven or calico material
- Calico bandage can be washed and reused, so is ideal for training
- Bandages are individually wrapped

Description	Size	Pack	Ref
Non - woven	90cm x 90cm x 127cm	1 x 50	9560
Calico - hemmed	90cm x 90cm x 127cm	1 x 50	9562



Ref 9560



Ref 9562

Steroswab

Sterile Non-woven Swabs - 4 ply
Help cushion, protect and prevent dryness in an open wound

- Used as an additional absorbent layer for bodily fluids
- Applied in conjunction with a low adherent dressing
- 4 ply cotton gauze swabs available in three sizes

Description	Pack	Ref
Steroswab - 5cm x 5cm	25 x 5	1903
Steroswab - 7.5cm x 7.5cm	25 x 5	1904
Steroswab - 10cm x 10cm	25 x 5	1905



Steroswab

Non-Sterile Gauze Swabs - 8 ply

- Used as an additional absorbent layer for bodily fluids
- Applied in conjunction with a low adherent dressing
- 8 ply cotton gauze swabs available in three sizes

Description	Pack	Ref
Steroswab - 5cm x 5cm	1 x 100	1900
Steroswab - 7.5cm x 7.5cm	1 x 100	1901
Steroswab - 10cm x 10cm	1 x 100	1902

