


Sterostrip-Sutures

Are the strongest, most reliable sutures on the market.



- The nylon strengthening fibres run the length of each suture adding to the strength of the suture, and preventing them from snapping under stress through normal body movements
- Prevent wounds from re-opening, ideal for knife cuts and similar deep wounds
- The medical grade adhesive is formulated to ensure the strip stays on longer, which in turn leads to a better healing wound or incision
- Hypoallergenic for less skin reactions
- Individually wrapped and sterile

Description	Size	Pack	Ref
Sterostrip - Sutures	3mm x 75mm	10 pouches of 5	7858
Sterostrip - Sutures	6mm x 75mm	10 pouches of 3	7859