



**1m tubular gauze bandage with applicator**

9540/1M



# Sterogauz Tubular Bandage

The ideal tubular bandage for holding dressings in place on the body.

- Ideally suited to apply over areas of the body being treated with creams and ointments, as it prevents transfer onto clothing and bedding
- Has natural stretch which provides comfort and conformability around body contours
- A range of sizes suitable for the whole body

Description	Size	Pack	Ref
Sterogauz 20m - dispensing carton	00 (1cm)	1 x 6	9545
Tubular gauze bandage	01 (1.5cm)	1 x 6	9540
	12 (2.5cm)	1 x 6	9541
	34 (3.5cm)	1	9542
	56 (5cm)	1	9543
	78 (8cm)	1	9544
Plastic applicator		25	9551/25
Plastic applicator		100	9551/100
<b>1m tubular gauze bandage with applicator</b>	01 (1.5cm)	1	9540/1M

# Roll-on Finger Bandage

- The finger bandage can be easily applied to either fingers or toes
- The unique rolled design removes the need for applicators
- Easy to use finger bandage and cots

Description	Colour	Pack	Ref
Finger bandage	White	1 x 6	3061WT
Finger bandage	Blue	1 x 6	3061BL



1. Put roll-on finger bandage on end of finger



2. Roll the bandage down finger



3. Then twist and roll down second layer of bandage



4. Finished application