

Trauma Dressing



“ No dressing is quicker to use, or more effective in operation ”



Remove from waterproof pack and sealed sterile bag.



Place directly onto wound.



Speed of application is easy, plus direct pressure applied.



Direct pressure quickly stems blood flow, and protecting wound from contamination.



Quick to secure with Velcro fastening.



Casualty can now be moved, allowing medics to deal with other injuries.

Trauma Dressing

Applies instant direct pressure – quickly stems major blood loss.

Multiple casualties can be treated quickly. No dressing is quicker to use, or more effective!



- Quickly stems major bleed injuries
- Velcro fastening ensures end of bandage is secured quickly
- Can be used on any part of the body
- Two size options 18 x 10 and 18 x 15cm
- Vacuum packed to condense bandage for storage
- Sterile to prevent infection

Super absorbent pad
Absorbs over
Twice its weight
in fluid

Description	Ref
Trauma bandage 10cm x 18cm pad	6830
Trauma bandage 15cm x 18cm pad	6831

☎ To order call **0161 902 3030**
 ✉ email **marketing@steroplast.co.uk**
 🌐 or visit **www.steroplast.co.uk**

