

S.T.A.T. Tourniquet



The S.T.A.T. stops bleeding fast (Smart Technical Application Tourniquet)

5 seconds to apply! Treat dozens of casualties in minutes!

Quick and easy to use on adults and children. A precise mechanism allows for 2mm adjustments.

No skin reactions – it's **latex-free**

The perfect tourniquet for mass casualty incidents. Any member of the public can help out in an emergency.

Each tourniquet has visible, clear instructions printed directly on the tourniquet.

The S.T.A.T. Tourniquet design is **similar to that of an everyday cable tie**, so it's simple and intuitive, and the reason it only takes 5 seconds to put on a casualty.

The highly visible **child-proof** release lever stops the tourniquet coming off, but also allows for ease of removal once the individual is receiving medical attention.

The **"Timestrip"** self-counting liquid timer helps triage victims. Allowing the individuals applying the tourniquet to time exactly how long it's been on.

Rounded edges are kinder to skin and soft tissue resulting in less trauma to the patient. The finger loop in the strap allows the user to get a strong grip during application, to apply the right amount of pressure.

Should the situation call for it, two S.T.A.T.'s can be combined and used alongside a trauma dressing to create a torso compression bandage, providing additional compression over a deep wound or laceration.

S.T.A.T. Tourniquet Orange - single use tourniquet - Code 1146

S.T.A.T. Tourniquet Blue - Colour coded for first response training - Code 1146/BLU

CE Marked



Fits in the palm of your hand



Self-counting liquid Timestrip



Child-proof release lever

Order your S.T.A.T. Tourniquet now, by quoting 1146 or 1146/BLU

To order call 0161 902 3030

email marketing@steroplast.co.uk

or visit www.steroplast.co.uk



Setting Standards Across Seven Sectors

Emergency

Hospital

Sport

Veterinary

Medical

First Aid

Retail