

Sterotracs Skin Traction Kit



Skin Traction Kit

Gain and maintain bone alignment after dislocation, to correct mild deformity and relieve pain.

Adhesive Skin Traction Kit

The extension strapping in the Sterotracs kit is coated with a hypo-allergenic adhesive, making it ideal for use on patients who react to zinc oxide or rubber based adhesives.

The light-weight retaining bandage is specially designed for the kit.

It comes with a soft foam strip which protects the foot, and a spreader plate with weight carrying cord for use either with fixed or balanced traction.

Non-adhesive Traction Kit

A specially developed ventilated traction bandage is used in the Sterotracs kit which, in conjunction with the high-quality retaining bandage, eliminates the need for adhesive extension strapping.

It comes with a soft foam strip to protect the foot, plus a spreader plate with weight carrying cord for use either with fixed or balanced traction.

What Is Sterotracs used for:

- Skin traction can be used to gain or maintain bony alignment.
- As in a dislocation or fracture.
- To secure immobilisation of an inflamed or injured joint.
- To correct mild deformity.
- To relieve pain.
- It is a valuable means of applying traction to a limb and is free from certain risks associated with skeletal traction.
- Skin traction may be used for long-term treatment or until such time that surgical procedures have been determined.

Description	Pack	Ref
Sterotracs - Adult	18	9517
Sterotracs - Child	24	9518
Sterotracs adhesive - Adult	18	9515
Sterotracs adhesive - Child	24	9516