

Healthcare Buyers Guide



Version 14B

2 NEW WEBSITES

Browse dozens of new products & numerous solutions to your sourcing problems!

When you're looking to save time & money, buying tried and tested **specialist trauma products**, **advanced life support** or a reliable comprehensive **range of medical essentials** our 25 years manufacturing experience is invaluable.



www.sterosport.co.uk



www.steroplast.co.uk



5 NEW CATALOGUES



Scan the QR code to view online or contact us to request your free copies.

NHS APPROVED PRODUCTS

100% SATISFACTION GUARANTEED

25 YEARS RETAIL EXPERTISE

Bags 50-53

Rucksacks

Waterstop® Specialist Bags

Bandages | Tubular 46-47

Sterogrip Bandages

Sterogauz Bandages

Finger Bandages

Bandages 45

Conforming Bandages

Crepe Bandages

WOW Bandages

Boxes, Stations & Cabinets 54-55

Defibrillators 18

Diagnostics 16-17

Disposables & First Aid Essentials 58-61

Miscellaneous items

Dressings 14

Ambulance Dressings

First Aid Dressings

Trauma Dressings 15



Public Access Trauma Kits 10-11
a.k.a "PACT Kits"



First Aid Kits | Full Range 4-5
Specialist First Aid Kits

First Response Kits 32-33

For "Grab-and-Go" Emergencies

40 & 70 Piece Kits

Hot & Cold Therapy 44

Hot & Cold Packs

Ice Packs

Manikins 19



Plasters & Wound Care 38-41
Includes Which? Recommended Plasters



Pharmaceuticals
Visit our website for all pharmaceutical products

Scissors

First Aid Scissors

Tweezers

57

View all catalogues on our homepage...

HSE at work Showcase

British Standard & HSE compliant First Aid Kits & supplies for common workplace injuries.



Scan the QR code to view online



Sterosport Buyers Guide

Sports tapes, strappings & first aid essentials specially selected for the world of sport.



Scan the QR code to view online



Retail Buyers Guide

A range of high-street healthcare products made with 25 years of manufacturing expertise.



Scan the QR code to view online



ESS Buyers Guide

A selection of the very best brands and products for the Emergency Services.



Scan the QR code to view online



Veterinary Buyers Guide

NHS quality dressings, wound care and specialist kits for the veterinary sector.



Scan the QR code to view online



Sterilising Fluids

& PPE

58-59

Tapes

Zinc Oxide Tape

Sports Tape

Sleek Tape

42-43

Training & Resus

Personal Protection Pouch

One Way Valves

20-21



Wipes

PDI Range


58-59

Contents Part 2 | See overleaf →

For our Emergency Services products, visit www.steroplast.co.uk

or contact us to request a **FREE** copy of the **Steroplast Emergency Services Brochure**.

Call: 0161 902 3030 | Email: info@steroplast.co.uk





Available in Small, Medium & Large

BS8599-1 Workplace First Aid Kit

24-25

Ensuring the well being of staff in the workplace. Keeping companies Health and Safety compliant.



Also available in a bag

Childcare First Aid Kit

30

Put together in consultation with the National Childminding Association, the Health & Safety Executive and Pre-school Learning Alliance.



Also available in a bag

Lone Worker/Off-Site First Aid Kit

28

Equipped with contents suitable for people working alone so that they can cope with a variety of accidents.



Includes S.T.A.T Tourniquet

Steroplast Bleed Control First Aid Kit

12-13

Supplied with specialist tourniquets for catastrophic bleeding.



Sizes 1-5, 1-10, 11-20 & 21-50 persons

HSE - Health & Safety First Aid Kit

26-27

The ideal kit for a low-hazard environment, compliant with Health & Safety law.



Great Value First Aid Kits

40 & 70 Piece First Aid Kits

33

Same high quality contents as all other first aid kits in a compact nylon case with see through compartments.

In an emergency, our First Aid Kits make finding products quick & easy...

We've got a range of kits specifically designed for your individual needs.



Eyecare 22-23

Emergency Eyecare Kits

Burncare 34-35

Emergency Burncare Kits

Catering 36-37

Specialist Catering Kits



Specialist
Kits

Biohazard & Sharps Disposal 48-49

The safest way to deal with sharp objects and infectious bodily fluids.



BS8599-2
Compliant

Vehicle First Aid Kits 29

Steroplast's British Standard vehicle first aid kits are supplied with Trauma bandages for life threatening wounds.

Taxi (PSV) First Aid Kits 29

Drivers of taxis/Public Service Vehicles (PSV) must have a first aid kit onboard in case a passenger gets injured.

You can rest assured that everything we supply is of the highest quality, is tried and tested, comes complete with robust guarantees, gets delivered on time, and has been designed to meet the toughest UK, EU and global standards.



Open the camera app on your smartphone device and hover over the QR code to find out more

"We know from what you've told us that you like the fact you can always talk to the same person, that knows you, and knows what you need for your business"

At Steroplast you'll always be greeted by a friendly voice!

No long-winded phone menus to navigate, when all you want is to talk to someone who is friendly, courteous and helpful!

Fast service every time

We won't keep you waiting. All our branded products are available ex-stock, so delivery is fast - within **24 hours**

in the UK, and within just a few days for Europe and the rest of the world.

Large back-up stocks for fast deliveries

We hold **large back-up stocks** to ensure faster fulfilment especially during seasonal demand surges e.g. sport-related products. In the unlikely event of a delay in delivery, you'll always be kept informed.

Ordering made simple...



0161 902 3030
+44 (0)161 902 3030



Ordering online...

www.steroplast.co.uk
www.sterosport.co.uk



Order via email:
sales@steroplast.co.uk



Steroplast Healthcare Ltd
FREEPOST NEA 16273
Manchester M22 3EA

Same day pick-up

If you're in the Manchester area you can collect your order the same day.

Next working day delivery

Orders tracked on-line: you'll know where your order is at all times.

Drop shipments: If at any time you would like us to deliver directly to your customer, we can do that. Sometimes people choose this for speed, or just because it's easier.

Payment method

Payments can be made by **bank transfer, credit card or debit card.**


Instant credit account available.

Friendly after sales care:

You'll find our after sales care is professional and courteous. The customer service team are always here to help.

From time to time something may not run quite as smoothly as we or you would like; that's when you'll find out how much care we take in providing a solution, to ensure your complete satisfaction.

Did you know Steroplast can supply all your pharmaceuticals and medicines? You can **now** buy all your pharmaceuticals and medicines from one location saving you time and money.



Prescription Only Medicines (POM)



Pharmacy Medicines (P lines)



General Sales List (GSL)

We can supply POMs to holders of a valid Wholesale Dealer's Licence, OR on receipt of a 'POM Authorisation Form.

Your options are:

- 1) A POM authorisation form signed by a Doctor
- 2) Through a registered nurse authorised by a Doctor
- 3) Through a registered paramedic

Steroplast Healthcare is licensed by the MHRA

To sell a comprehensive range of ***POM, *P lines and *GSL** products.



POM - Prescription Only Medicines

Any medicine classified as a **POM is restricted** for use under medical supervision on a prescription only.

P lines - Pharmacy medicines

Medicines classified as **"P" lines are restricted** for use under the supervision of a pharmacist.

GSL - General Sales List products

Available for sale in general retail outlets (**restricted by pack size**). These items are considered safe for use without medical supervision.

With excellent prices on offer, and an opportunity to save time and money, plus a wide range of new, unique and innovative products, calling Steroplast could be the best thing you do to-day! **T: 0161 902 3030 E: sales@steroplast.co.uk**

Through our global manufacturing partnerships developed over 25 years, we can offer you new, unique and innovative products designed for 7 individual markets, **delivered to your door within 24 hours.**

25 years
Manufacturing
Experience &
Expertise

Having global suppliers guarantees you get great prices!

We have gone to the ends of the earth, literally, to ensure we can deliver major price advantages to you.



100% satisfaction guaranteed: Zero risk to you.

We're so confident in the products we offer, we guarantee every purchase you make. If for any reason you're not entirely satisfied with any part of your order, once delivered, just let us know and we'll put it right, no matter what it is. Starting a partnership with Steroplast is **RISK FREE!**

Working in partnership, helping your business grow

We recognise the importance of a true partnership in business, and have come to rely on partnerships with suppliers around the world. We're here to offer help and guidance when required.

Sourcing new and bespoke medical product solutions for the **British Red Cross**

Steroplast have been working with the **British Red Cross** for a number of years as their exclusive supplier of medical products. We have also helped them overcome many of their logistical challenges in sourcing new and bespoke products tailored to their individual requirements.

We can do the same for you.

Take this opportunity to contact us and explore ways we can work together to provide you with a solution to your medical product sourcing problem.

Many of the specialist products we supply to the **British Red Cross** and throughout the emergency services can be found in our **"Emergency Services Sector"** brochure covering the following categories:

- Airway Management
- Diagnostics
- Evacuation/Immobilisation
- Injection/Infusion
- Pharmacy
- Trauma

Request your **FREE** Emergency Services Catalogue

For more **detailed information** go to our new **Homepage** and search under **Browse Categories.**

Click on Sectors for
Emergency Services
Equipment

OR

Advanced Life
Support



PLUS you can take advantage of your first time buyers special offers and discounts when you make contact. Call **0161 902 3030** or email **info@steroplast.co.uk**

From the UK, to 57 countries world-wide.

Exports account for more than 25% of our turnover, demonstrating our competitive advantage.



"I was very impressed with all the customer service staff that I dealt with. A seamless service. Thank you"

IH Physio Solutions

"Fast efficient service, paid for the order and was here the next day. Very friendly and polite staff, excellent customer service"

Simmtronic

"Excellent communication and customer service"

Ignition Search

"You're doing a great job during difficult times - thanks!"

Indigo

"Great company to deal with and go out of their way to support us"

SMS Pharmacy

"Fantastic customer service and availability always. Would highly recommend"

Barrow AFC

An ethical business offering unrivalled quality

At Steroplast we are committed to the very **highest standards of quality and ethical business practice**. We abide by the most stringent medical and ethical guidelines monitored by each of the following governing bodies:

Supplier Ethical Data Exchange



We are members of Sedex an organisation for businesses committed to continuous improvement of their ethical supply chains.

The Medical Device Directive



Steroplast conducts its operations according to the requirements of **BS EN ISO 13485:2012** supported by:

The Medical Device Directive (M.D.D.) 93/42/EEC and The Wholesale Dealers Authorisation Licence.



Our quality system is fully compliant with the **Medicines and Healthcare products Regulatory Agency**.

Give us the chance to prove the claims we make, and help you succeed. Get in touch today and find out how we can help. Email us on **info@steroplast.co.uk**

If a picture's worth a 1000 words, then a video must be worth 10,000!

Learn more about us and our products in high quality, engaging short video clips: 30 – 90 seconds each.



For more information about Steroplast and our product range, have a look at our videos on our home page.

Public Access Trauma (PACT) Kit



What is a Public Access Trauma Kit (PACT)?

The Public Access Trauma Kit was designed for non-medically trained people to support life-threatening injuries such as a major bleed or someone who is not breathing.

In response to the devastating Manchester Arena bombing in 2017, the National Counter Terrorism Security Office collaborated with medical experts to create this kit, which has now become the government standard.

Advantages, features and benefits...

- Trauma Dressings designed on Military style dressings for major bleeds
- For use by the general public, including young people and people with no previous medical or first aid training
- Products allow for self-treatment of the walking wounded
- To help businesses, public places etc be prepared for mass casualty incidents
- The contents can treat major bleeding, unresponsiveness and irregular or absent breathing

What is in a PACT Kit?

A basic kit, contains 2 windlass tourniquets, 4 trauma dressings, 2 face shields for CPR, 2 triangular bandages, 4 PPE packs and a simple to follow instruction guide.

Who should have this kit?

All public accessible locations (PALs) should have one or more kits. This could be a stadium, arena, museum, bar, restaurant etc.



Public Access Trauma (PACT) Kit

DOES A PACT KIT REPLACE A FIRST AID KIT?

NO, all businesses with over 5 employees must still meet their legal obligations for workplace first aid.



How many people does a PACT Kit treat?

A standard kit can treat 2-4 people depending on the severity of the injury. Talk with us about larger kit options and enhanced kits.

How many do I need?

This will depend on your size and risk assessments whilst also considering your vulnerability and counter terrorism plan.

Suggested government guide below)...



SUGGESTED LOCATIONS FOR PACT KITS

(Source: ProtectUK.Police.UK)

Business Type	Main Entrance	Fire Exit	Concourse	Muster Point	Refuge Area	Ticket/Booking Hall	Platforms	Drivers Cockpit
Offices	✓	✓	✓	✓	✓			
Retailers	✓	✓						
Stadiums/Arenas/Hotels/Conference Centres/Museums	✓	✓	✓	✓	✓	✓		
Transport Locations/Vehicles	✓	✓	✓			✓	✓	✓
Public Areas Including Parks/Play Areas	✓			✓				
Events Including Fates/Festivals/Shows	✓			✓		✓		

Bleed Control Kit

Only 5 seconds
to apply



The S.T.A.T. Tourniquet stops bleeding fast. Supplied with Sterostop Bleed Control kits.

The perfect tourniquet for mass casualty incidents. Any member of the public can help out in an emergency.



3 min

A casualty who is bleeding from an artery can die in as little as **3 minutes**

5 seconds to apply! Treat dozens of casualties in minutes!

Quick and easy to use on adults and children. A precise mechanism allows for 2mm adjustments.

Design is **similar to that of an everyday cable tie**, so it's simple and intuitive, and the reason it only takes 5 seconds to put on a casualty.

The highly visible **child-proof** release lever stops the tourniquet coming off, but also allows for **ease of removal** once the individual is receiving medical attention.

Unique "Timestrip" self-counting liquid timer helps triage victims. Allowing the individuals applying the tourniquet to time exactly how long it's been on.

Should the situation call for it, two S.T.A.T.'s can be combined and used alongside a **trauma dressing** to create a torso compression bandage, providing additional compression over a deep wound or laceration.

S.T.A.T. Tourniquet Orange - single use tourniquet (Smart Technical Application Tourniquet) Order code 1146.

S.T.A.T. Tourniquet Blue - For first response training. Order code 1146/BLU.



Ref 6830

Trauma dressing - applies instant direct pressure – quickly stems major blood loss. Quick to secure with Velcro fastening.



Ref 6824

Sterostop - bleed control kit provides everything needed to stem the flow of blood, in an emergency.



"Timer starts counting once activated"
Rounded edges avoid skin cuts, less trauma to the patient.



No skin reactions
it's **latex-free**.

Ref 1146

Each tourniquet has clear instructions printed on the tourniquet.

Critical Injury Kit

Made to British Standard BS8599-1

It's a matter of life and death



- Chito-Sam 100 Haemostatic Dressing x 2
- Foil Blanket
- First Aid Guidance Leaflet
- Gloves x 2 pair
- Trauma Dressing x 2
- S.T.A.T. Tourniquet

The perfect tourniquet for mass casualty incidents. Any member of the public can help out in an emergency.

- Tuff Cut Scissors 7"

Description	Pack	Ref
Critical injury kit	1	8599-C



Steroplast Critical Injury Kit

Specifically designed to control and stop severe blood loss in the event of a stabbing or similar serious injury.

The kit is equipped with tried and tested essential products that will enable you to save a friend, a colleague or loved one's life.

A Life saver designed for anyone to use.

The Steroplast Critical Injury Kit is supplied with easy-to-follow instructions so that anyone can help in a crisis, with or without medical training. It is specifically designed to control and stop severe blood loss in the event of a stabbing or similar serious injury, with everything you need to treat an injured victim and save a life. This is an ideal first aid kit for chainsaw users and construction workers working with high-risk tools and machinery.

Brightly coloured bag or box

Easy to find and identify in an emergency.

Trauma Dressing

Applies with instant direct pressure – quickly stems major blood loss.

Multiple casualties can be treated quickly. No dressing is quicker to use, or more effective!

- Quickly stems major bleed injuries
- Velcro fastening ensures end of bandage is secured quickly
- Two size options 18 x 10cm and 18 x 15cm
- Vacuum packed to condense bandage for storage
- Sterile to prevent infection



Absorbs over 7 times its weight in fluid

Description		Ref
Trauma bandage	10cm x 18cm pad	6830
Trauma bandage	15cm x 18cm pad	6831

Ambulance Dressings

Specifically for severe wounds to stem heavy blood loss.

Ambulance dressings are made with a thicker wound pad for greater absorbency.

Life saving dressings following serious trauma



Also available in 12's



Also available in 12's

Description	Size	Pack	Ref
No 1	10 X 12cm	1 x 72	7001B
No 2	20 X 15cm	1 x 48	7002B
No 3	20 X 28cm	1 x 30	7003B
No 4	20 X 32cm	1 x 24	7004B
No 7	Small - 5 X 6cm Finger	1 x 140	7000B
No 8	Medium - 10 X 8cm	1 x 72	4702B
No 9	Large - 13 X 11cm	1 x 48	4703B
No 16	Eyepad	1 x 156	1912B

Steropax First Aid Dressing

High quality wound dressings for inclusion in first aid kits.

Available with a paper wrapper



Also available in 12's



FINGER

Also available in 12's

Description	Size	Pack	Ref
LARGE	18 x 18cm (HSE spec)	1 x 48	1910B
MEDIUM	12 x 12cm (HSE spec)	1 x 72	1911B
EYEPAD	Boxed	1 x 165	1912BX
FINGER	3.5 x 3.5cm	1 x 100	1915B

Trauma

Whatever the type and size of the injury, we've got the ideal product to treat it quickly and effectively in an emergency.

- Bleed Control Kits
- Tourniquets
- Trauma Dressings
- Haemostatic Dressings
- Chest Seals



Talk to us about our Burnshield range
For 1st, 2nd and 3rd degree burns and scalds, dressings or gel start cooling burns instantly.



Prometheus ChitoGauze XR Pro
Controls moderate to severe haemorrhage wounds.



STAT Tourniquet

Simple to use, ideal for mass casualty, non medical trained use, with 60 minute timer.



SOF Tactical Tourniquet - Wide

Light, compact, robust, the medical professions choice.



Tactical Blast Bandage

For traumatic amputations, burns and large pattern wounds. Sterile.



Trauma Dressing

Quickly stems major bleeds, sterile, vacuum packed. Available in 2 sizes.



Celox

Offering the range of Celox haemostatic products to treat severe bleeding.



Corius Seal

Non sterile, adhesive hydrogel, covers and protects acute wounds, absorbs exudate.



Chito Sam

A high performance haemostatic product to treat major bleed. Simple and easy to use.



Russell Chest Seal

Fast, easy to apply to cover chest wounds.



Fox Seal - Vented Chest Seal

512 venting pathways, providing more ways to release pressure.

FREE

Request your FREE Emergency Services Catalogue

Steroplast has the very best brands and products for the Emergency Services under one roof. With a comprehensive range of medical products that service your needs are all the benefits of working from one location.

Injection/Infusion

Any type of medical injection or infusion.

- Insulin
- Vaccines
- Blood products
- Chemotherapy
- Anticoagulants

Diagnostics

Our high quality diagnostic equipment includes:

- X-ray
- Ultrasound
- ECG
- Blood pressure monitors
- Stethoscopes
- Thermometers
- Pulse oximeters
- Blood glucose monitors

Pharmacy

Our pharmacy is staffed by qualified pharmacists and offers a wide range of pharmaceuticals. We also offer a range of medical supplies and equipment.

Emergency Services Sector
Better quality, more choice, one location

Diagnostics

Here is a selection of the diagnostic equipment available at Steroplast, visit our website or see QR code for full list.



Open the camera app on your smartphone device and hover over the QR code to find out more

Forehead thermometer

Available in 25's and 50's

The unbreakable and easy to use forehead thermometer gives fast and accurate readings.

Advantages, features & benefits

- Completely safe
- The pocket-sized product can be taken anywhere
- Fast and simple to use
- Temperature range: 35 - 40 Celsius



Ref 1811/25
Ref 1811/50

- Thermometer dimensions: 89mm x 19mm
- Latex-Free
- Hospital and home use
- The fast and simple temperature monitoring system can be used as pre and post-operation surveillance.

Digital ear thermometer

Compact, ergonomic digital thermometer with a 1 second measurement response time.

- Holds a capacity of 10 memory records
- Auto shut off
- Approx. 4,000 uses between batteries
- Fast and accurate



Ref 1128



Ref 1802

Braun ThermoScan ear thermometers

Ear thermometer measures the infrared heat generated by the eardrum and surrounding tissue.

- Dual scale (celsius or fahrenheit)
- Ultra soft and gentle tip is warmed for extra comfort



Ref 1817



Ref 1799

Digital thermometer

A safe, accurate and simple way to check and monitor your body temperature.

- Accurate and easy to read display
- Simple to operate
- No glass or mercury
- Provides fast and accurate body temperature reading
- Large screen for fast and easy reading of results
- Automatic shut-off function preserves battery life



Ref 1818

Berrcom non-contact infrared thermometer

This forehead thermometer is designed to take a temperature with a 3-5cm distance from the forehead. It uses infrared technology to measure temperature and is quick and easy to use.

- Non-contact for hygienic, non-intrusive use
- Audible fever alarm will alert you to a temperature of over 38.6°C
- Easy to read, backlit LCD screen
- Temperature displayed in either Fahrenheit or Celsius
- Body measurement range: 25°C - 42.9°C
- Auto power off after 30 seconds



Ref 1796

Twin tube Sprague Rappaport stethoscope

Traditional design stethoscope with two tubes connecting the chest piece to the head frame.


- Two tubes connecting the chest piece to the head frame
- Supplied with interchangeable chest pieces for adult and paediatric use
- Available in 6 additional colours Red/Purple/Pink/Orange/Blue/Green/Black

Ref WSTE0116/RED



Cardiology stethoscope

- Stainless steel double head chest piece with a deep cone-shaped bell gives excellent acoustic sensitivity
- The thick tube offers resistance to ambient interference
- Binaural with reinforced stainless steel inner-spring for durability and portability
- 3-year warranty



Ref 1063

Deluxe cardiology stethoscope

- Excellent sound conduction
- Zinc alloy chest piece
- Flexible binaural for optimal user comfort
- Durable soft silicone ear tips
- Supplied with spare ear tips
- 3-year warranty



Ref 1118

More Stethoscopes available on our website www.steroplast.co.uk

Doctors surgeries and hospitals within a 20 mile radius from our office, can have products delivered within 2-3 hours



Ref WSTE0077

Creative medical fingertip pulse oximeter

The tough, reliable and accurate PC-60B1 EMS Professional Model.

- 3 Year drop warranty - Super tough EMS/Military model
- Automatic ON/OFF for ease of use
- 4 Direction dual colour LCD display
- Plethysmographic waveform allows optimum finger placement
- Large backlit display of SpO2, PR, PI and pulse bar
- Supplied with carry pouch, lanyard and batteries



Omron M6 AC LED blood pressure monitor

Clinically validated monitor incorporates the most advanced Omron sensors and indicators for accurate and reliable blood pressure readings.

- High performance, quality digital blood pressure monitor
- 100 reading, 2 user memory
- Body movement detector
- Hypertension tracker
- Clinical and Diabetic validation
- Stroke prevention



Monitor Cuff - Large (Ref WSTE0111)

Freestyle Optium Neo blood glucose meter

Designed to make life simpler for people who use insulin.

- Accurate blood glucose and blood ketone testing
- Insulin logging helps track insulin doses
- Fast, five second test time
- Notification when results need attention
- Clear and sharp screen
- Matte screen reads like paper, no glare
- Icon driven display

Ref FFR2029



Defibrillators

What is the difference between the semi-automatic and fully automatic models?

There is only one difference between the two AEDs. If in an emergency situation a shock was necessary, the semi-automatic model would prompt the user to press the button to deliver the shock.

Whereas, the fully automatic AED is designed to skip this step and automatically deliver the shock required.



iPAD-SP1 - Semi-Automatic Defibrillator
Solid and reliable defibrillator (Ref 1844)

iPAD-SP1 - Fully-Automatic Defibrillator
Simple and fully-automatic (Ref 1875)

Zoll AED Plus Lay Responder Defibrillator
Only AED with 'Real CPR Help' (Ref ZOLL0001)

Schiller FRED Easyport AED
Small and lightweight semi-automatic (Ref 1729)



iPAD SP1 - Child and Adult Replacement Pads
For both adults and children (Ref 1845)

iPAD SP1 - Replacement Battery
For fully-automatic and semi-automatic iPADs (Ref 1855)

Zoll AED Plus - Adult Stat-Pad II
Zoll AED plus defibrillator pads (Ref ZOLL007)

Schiller FRED Easyport AED - Adult Pads
Pre-gelled in a sealed pack (Ref 1733)



iPAD-SP1 - Defibrillator Training Unit
Functions an instructor needs (Ref 1849)

Defibsafe 2 External Lockable Cabinet
Superior protection for defib kept outside (Ref 1882)

Zoll AED Plus - Trainer 2
Features and functions an instructor needs to teach students (Ref ZOLL0002)

Indoor Defibrillator Cabinet with Strobe Light and Alarm
Secure indoor storage (Ref 1876)

See website for defibrillators accessories and other products www.steroplast.co.uk

Manikins

We supply a range of Manikins – **Laerdal**, **Brayden**, **Ambu** and more.

The perfect first aid tool to create an ideal training experience.

This method of physical hands on learning is necessary for medical training as it helps students to understand the procedure first hand, allowing them to practice techniques on a life-like torso which simulates patients.



Brayden Pro CPR Manikin

For data-driven performance-focused debriefing of CPR training. Brayden Pro & Tablet (Ref BRAY25787)



Brayden CPR Manikin

This anatomically correct Brayden CPR training manikin (Ref BRAY25785)



Choking Charlie

Training tool for teaching the Abdominal Thrust Manoeuvre (Ref 1350)



Brayden CPR Training Manikin with Red Indicator Lights

Manikin (Ref BRAY25780)



Laerdal Little Baby - Q CPR Training Manikins

Laerdal Little Baby Q CPR Light Skin Manikin (Ref LEAR010)



Laerdal Little Anne - Q CPR Training Manikins

Laerdal Little Anne Q CPR Light Skin Manikin (Ref 1094)



Laerdal Junior - Q CPR Training Manikins

Laerdal Junior - Light Skin Q CPR Training Manikin (Ref 1096)



Laerdal Little Baby - Q CPR Training Manikins

Laerdal Little Baby Q CPR Dark Skin Manikin (Ref 1106)



Laerdal Little Anne - Q CPR Training Manikins

Laerdal Little Anne Q CPR Dark Skin Manikin - 4 Pack (Ref L3050)



Laerdal Junior - Q CPR Training Manikins

Laerdal Junior - Dark Skin Q CPR Training Manikin (Ref 1108)

CPR

Pocket Mask

The mask that makes resuscitation safe for the provider, the one-way valve with membrane prevents cross infection.

- Supplied in clamshell case, portable and convenient
- Improves the effectiveness of cardiopulmonary resuscitation

Description	Pack	Ref
CPR pocket mask	1	9598



Resusciate

Face Shield

The safest way to give mouth to mouth resuscitation.

Offering protection to both rescuer and casualty.

Description	Pack	Ref
Individually packed resuscitation face shield (see page ????????)	1 x 100	9589
Resuscitation face shield	Roll x 36	9590



Available on a roll of 36

- Used by CPR trained rescuers
- For inclusion with corporate CPR training courses, and in first aid kits
- Includes printed on CPR instructions



Training Packs

(Candidate Packs)

A high quality first aid training pack as used by numerous training companies, contents of which are:

- Medium training bandage
- A pair of gloves
- Face shield,
- Plaster
- Hemmed triangular bandage
- Safety pin



Description	Pack	Ref
First aid training packs	1	8698
HSE first aid at work guidance leaflet	1 x 20	LEAFIR007
Safety pins 6's assorted	1 x 100	1518/100

Training Dressing

Non-sterile extra thick training bandage



Description	Pack	Ref
Traning bandage - 12cm x 12cm	1 x 15	FFK1001

Resusciade

One Way Valve

The mouth to mouth resus that allows hygienic contact with a patient whilst resuscitating them.

- Used and promoted by CPR trained rescuers
- Ideal for corporate CPR training courses and for first aid kits
- Offers protection for both rescuer and casualty when used in CPR
- Includes printed on CPR instructions
- One way valve
- Single use



Resus - one way valve

Keyring Face Shield

- Available with your own logo and branding
- **Contact us for more information**



9583/25

Description	Pack	Ref
Resus – one way valve	1 x 25	9582/25
Resus – with ear loops	1 x 25	9583/25
Keyring face shield	1 x 25	9567/25

Personal Protection Pouch

A compact personal protection kit for first aiders to give CPR without risk of infection.


Easy to store compact kit, ideal for first aiders

- Can fit easily in your pocket or bag to be used at anytime
- Velcro fastening, for quick and easy access
- Made from strong nylon material
- Can be attached to a belt for hands free access



Contents: apron, gloves, resus and wipes

Description	Pack	Ref
Personal protection pouch	1	8132



Eyewash Station

- Eyewash should be wall mounted where there are toxic, corrosive, acid, or caustic chemicals
- Perfect to remove sand and dust particles which can damage the eye if not removed quickly

500(W)x750(H)

Ref 1422

The first 10 seconds are critical


The first 10 seconds after the eye has been exposed to a hazardous substance are critical, especially when that substance is corrosive.

Eyewash stations must be fixed to a wall near to the potential hazard. e.g. next to a fork lift charging point as battery acid splashes could cause blindness. Eyewash should be wall mounted where there are toxic, corrosive, acidic, or caustic, chemicals.

Sterowash is perfect to remove sand and dust particles which can damage the eye if not removed quickly.





Ref 1200




EMERGENCY EYEWASH STATION

& guidance for use







STEP 1: Twist off the top of the pod or bottle to open the eyewash solution provided.




STEP 2: Open eyes wide and position the open bottle of solution above the eye. Use the mirror provided if using without assistance.




STEP 3: Gently squeeze the bottle allowing the liquid to flow freely into and across the eyeball.



STEP 4: Use the mirror provided to observe the eye. Seek further medical assistance if necessary and use eyepad dressings to protect the eye.



STEP 5: Discard of the bottle after use. DO NOT REUSE. Replace any used contents of the kit for use in the next emergency.



Once fixed on the wall where it'll be seen every day, it's there permanently in case of emergency.

**Sterile
Eye Wash
& Wound
Wash**

**Now
On Drug
Tariff**



Sterowash (20ml)
1 x 25 pods

Ref 1400



Eye & wound wash station
(20ml x 10)

Ref 1410

Sterowash saline 20ml pods treat both eyes and wounds and are available on Drug Tariff, e.g. as prescribed by your doctor.



Eye wash station (500ml x 2)

Ref 1416

The 2-in-1 eye and wound wash could save someone's sight!

In an emergency a twist off top makes Sterowash ready for use in a second!

It's also perfect to cleanse wounds to reduce the risk of infection, as its sterile.

With a slimline neck for improved directional accuracy, and a twist off top ready in the blink of an eye, these pods are the only choice in an emergency.

**Now
On Drug
Tariff**

Drug Tariff

PIP Code	Description	Item	Ref
348-9168	Sterowash – sterile saline 20ml pods	10 x 25	1400W
348-9168	Sterowash – sterile saline 20ml pods	1 x 25	1400W/1

Class IIA Medical Device



Ref 9115

Sterowash Eye & Woundwash Kit

with wall bracket and 2 x 500ml bottles & 5 x 20ml pods.

The Ultra Durable Premier Box

Ideal for warehouses, factories, garages etc, as it protects the eye and wound wash, and keeps it clean. These stations are specifically designed for your workplace, you decide which kit you prefer.

Description	Item	Ref
Sterowash Eyecare Kit 2 (2 x 500ml)	1	9114
Sterowash Eyecare Kit 3 (3 x 500ml)	1	9116
Sterowash 500ml bottles - Shelf life 4 years 112 cases per pallet	1x10	1405
Sterowash 20ml pods - Shelf life 3 years 36 cases per pallet	10x25 1x25	1400W 1400W/1



Ref 9116

Sterowash Eyewash Kit 3

with (3 x 500ml) bottles of eyewash with eyepads provided to cover the eye while waiting for medical attention.

Workplace First Aid Kit

Compliant to British Standard BS8599-1



Open the camera app on your smartphone device and hover over the QR code to find out more

3 sizes available



Workplace First Aid Kits

Contents	S	M	L
Bandage - conforming	1	2	2
Burn dressing	1	2	2
First aid emergency leaflet	1	1	1
Foil blanket	1	2	3
Gloves - nitrile in pairs	6	9	12
Micropore tape	1	1	1
Plasters w/p	40	60	100
Resus device/face shield	1	1	2
Safety Pins	2	4	6
Scissors - tuff-cut	1	1	1
Triangular bandage	2	3	4
Wall bracket	1	1	1
Wipes - saline for cleaning	20	30	40
Wound dressing - eyepad	2	3	4
Wound dressing - finger	2	3	4
Wound dressing - large	1	2	2
Wound dressing - medium	4	6	8

We've included a burn dressing which ensures you can treat a skin burn, reduce the pain without the need to keep the burn under running water.

Ultra Strong Bracket



Security Seal Ref 1221RE/50



Ref 9661

If your staff are split over more than one floor, you must have at least one kit per floor. First aid kits should be **positioned where people see them every day, so that in the event of an emergency, people will remember where to find them.**

In an emergency our NEW British Standard - Workplace First Aid Kits make finding products quick and easy.

The very latest BS8599-1 first aid kits are designed to be user friendly in an emergency.



As injuries at work can be life threatening you should be able to reach a kit within 60 seconds.

If your staff are split over more than one floor, you must have at least one kit per floor. First aid kits should be positioned where people see them every day, so that in the event of an emergency, people will remember where to find them.

First aid kits save lives and are required by law!

In the event of sudden injury, failure to provide sufficient first aid could have extremely serious consequences. Employers and managers should therefore ensure that an employee who is injured or taken ill at work receives immediate attention.

Choosing the right first aid kit is now as easy as 1, 2, 3.

See the selection chart below for guidance.

Decide if your work environment is:

- 1 **LOW HAZARD** or **HIGH HAZARD**.
- 2 **How many employees** there are in your workplace.
- 3 **Which of the three sizes** of kit you will need, e.g. small, medium or large. Place an order, put in a prominent position, done. It's that simple!

You can rest assured that everything we supply is of the highest quality, is tried and tested, comes complete with robust guarantees, gets delivered on time, and has been designed to meet the toughest UK and EU standards.

"First aid kit compliance made nice and easy by Steroplast Healthcare"

BS8599-1 Workplace First Aid Kit in an ergonomic Premier Box

We recommend having a small workplace kit in your sports clubs / facilities

Description	Pack	Reference
Small workplace first aid kit	1	8599-PS19
Small Refill	1	8599-RS19
Medium workplace first aid kit	1	8599-PM19
Medium Refill	1	8599-RM19
Large workplace first aid kit	1	8599-PL19
Large Refill	1	8599-RL19

Choosing the right first aid kit is now as easy as 1, 2, 3.

SELECTION CHART

Number of significant areas in your workplace	Degree of hazard	Number of employees	Minimum requirements
<p>A significant area is an individual floor or building. You should be able to reach a casualty within 60 seconds of an initial call for assistance.</p> <p>e.g. An office and a warehouse on the same site constitutes as 2 significant areas</p> <p>Each area needs to be assessed on its own merits.</p> <p>Sports clubs / facilities: we recommend 1 small kit on each floor</p>	<p>1 LOW HAZARD - Office (e.g. shops, offices and all sports facilities).</p>	2 Less than 25	3 At least 1 Small kit
		25 - 50	At least 1 Medium kit
		More than 50	At least 1 Large kit per 100 employees
	<p>1 HIGH HAZARD - Warehouse (e.g. light engineering and assembly work, food processing, warehousing, dangerous machinery, construction, harmful chemicals etc).</p>	2 Less than 5	3 At least 1 Small kit
		5 - 25	At least 1 Medium kit
		More than 25	At least 1 Large kit per 25 employees

HSE 1-10, 11-20 Person First Aid Kits

As injuries at work can be life threatening you should be able to reach a kit within 60 seconds.

If your staff are split over more than one floor, you must have at least one kit per floor.

First aid kits should be **positioned where people see them every day**, so that in the event of an emergency, people will remember where to find them.

We've made it simple, because it really doesn't need to be complicated.

We give you less choice, but **guarantee** that with a Steroplast kit you'll be compliant, and on the **right side of the law!**




Description

Ref

HSE 1-10 first aid kit
26.5cm x 27.5cm x 9cm

8018



Durable wall bracket
Ref 8119WB

These (Health & Safety Executive) HSE compliant kits are ideal for an office or low-hazard environment and comply with Health & Safety law.

Deciding which first aid kit you need is actually very simple. So many companies offer so much choice it becomes confusing.

There's a 1-10 person, and 11-20 person kit, which reflects how many people they're designed for. If your workplace has more than 20 people, then adding an extra 1-10 person kit will ensure you're covered for 30 people, and compliant with Health & Safety law.

Contents HSE kit	8018 1-10	8019 11-20
First aid emergency leaflet	1	1
Gloves - in pairs	1	2
Plasters w/p	20	40
Safety pins - pack	6	12
Triangular bandage	4	6
Wipes - alcohol free	10	10
Wound dressing - eyepad	2	4
Wound dressing - medium	6	9
Wound dressing - large	2	3



1 Person bum bag Ref 8130



1 Person Ref 8124



1-5 Person Ref 8129

We guarantee that with an HSE Steroplast kit you'll be compliant, and on the right side of the law!



1-10 Person in clear bag Ref 8015



1-10 Person Ref 8013

Description	Pack	Ref
1 Person first aid bum bag	1	8130
Refill	1	8124RP
1 Person first aid bag	1	8124
Refill	1	8124RP
1-5 Person first aid kit bag	1	8129
Refill	1	8129RP
1-10 Person first aid kit bag	1	8013
Refill	1	8018RP

HSE 1 Person, 1-5, and up to 1-10 people

1 Person first aid bum bag (8130)

Specifically designed for an individual. The 1 Person first aid kit is supplied in a green nylon bum bag. Worn around the waist, leaving your hands free to dress a wound.

The quality of the contents in these kits are exactly the same as those used in the BS8599-1 (British Standard) kits.

1-5 Person first aid kit bag (8129)

These kits are the ideal solution to look after up to five people.

1-10 Person first aid kit bag (8013)


These kits are the ideal solution to look after up to ten people. The advantage with this kit is it's supplied in a soft nylon bag, as opposed to a hard plastic container. There is no compromise on the contents of this kit.

Contents	1 person 8130	1-5 person 8129	1-10 person 8013
Dressing - large	1	1	2
Dressing - medium	0	3	6
Dressing - eyepad	0	1	2
Bandage - triangular	2	2	4
Gloves - pair (medium)	2	1	1
Safety pins	1 bag	1 bag	1 bag
Wipes - alcohol free	2	4	10
Plasters hypo-allergenic	10	10	20
Guidance leaflet	1	1	1

Lone Worker/ Off-Site Kits

Compliant to British Standard BS8599-1

The contents of the BS8599-1 Lone Worker /Off-Site kits (AKA Personal Issue first aid kits) have been carefully selected to enable the user to deal with a wide variety of accidents while 'lone-working' or 'off-site'.



"Steroplast's BS8599 kits are the best choice if you want to protect your employees, and comply with Health & Safety Regulations"

See Bleed Control Kit & Critical Injury Kit (P12 - P13)

Critical Injury kits contain three items specifically for severe blood loss caused by serious injuries:

- Tourniquet
- Haemostatic wound dressing
- Trauma dressing

They are meant for employees that are exposed to greater risk of serious injury.

For people working in high-risk environments first aid is essential to prevent loss of life. If you work with any of the following:

- Power tools
- Agricultural equipment
- Forestry related equipment
- Equipment associated with construction sites

Where there is a risk that someone could inflict serious wounds on themselves they must have access to "Critical Injury Kits" or "Bleed Control Kits"



Lone worker in a pouch

Ref 8127 1



Off-Site Plus lone worker in a bag
Contain 500ml eye & wound wash

Ref 8599-T 3

Description	Pack	Reference
Lone worker in a pouch	1	8127 1
Lone worker in a bumbag	1	8130
BS8599-1 first aid offsite kit - box	1	8599-T/BX 2
BS8599-1 first aid offsite kit - bag	1	8599-T 3
Refill	1	8599-T/RP

Off-Site Plus Kit

Contains 500ml of sterile water which ensures the injury doesn't become infected by using potentially contaminated water sources. It also has a burn dressing, which requires no water to cool a burn.



Off-Site Plus kit in a box

Ref 8599-T/BX 2



Choosing the right vehicle first aid kit.

The size of the first aid kits required for motor vehicles is determined by the category of hazard and by the number of potential passengers.

Number of Passengers	Vehicle Examples	Number and Size of First Aid Kits Required
1 to 3	Mopeds, motorcycles, motor tricycles and quadricycles	1 x small
1 to 8	Cars, taxis, commercial vehicles	1 x medium
1 to 16	Mini buses, small buses	1 x large
17+	Buses, coaches	2 x large*

*Multi-level vehicles require at least two kits

British Standard Vehicle First Aid Kits for UK motorists.

When it's a matter of life and death, you must be prepared.

- If an accident happens, you'll want to be prepared for potentially life threatening injuries
- Major trauma injuries require specialist 'Trauma Dressings' in order that heavy blood loss can be contained

Easy to store - and comes in 3 sizes, boxed or bagged.

The Vehicle First Aid Kit takes up very little room in your car and could be placed under the seat or in the boot, out of harms way.



You can order kits individually

Description	Size	Reference
Vehicle kit	Small	8599MV-S
Vehicle kit	Medium	8599MV-M
Vehicle kit - Bag	Medium	8599MVB-M
Vehicle kit - Bag	Large	8599MVB-L
Taxi / PSV kit	Small	8160



BS8599-2:2019 Vehicle Kits

Contents	S	M	L
Burn dressing	1	2	4
First aid emergency leaflet	1	1	1
Foil blanket	0	1	3
Gloves - nitrile in pairs	1	2	5
Plasters w/p	5	10	20
Resus device/face shield	1	1	2
Scissors - tuff-cut	1	1	1
Triangular bandage	0	1	2
Wipes - saline for cleaning	5	10	20
Wound dressing - adhesive	0	1	2
Wound dressing - trauma L	0	0	1
Wound dressing - trauma M	1	1	2
Wound dressing - first aid M	0	1	2

Childcare First Aid Kit

You should be able to get to an injured child within 60 seconds.

Schools must have multiple kits on each floor.

Childminders, Schools and Nurseries

/// Contains hypo-allergenic plasters which are ideal for use on children ///

Quick-grab classroom first aid kit
Ref 8196

40 piece first aid kit
Ref 8133



70 piece first aid kit
Ref 8128



Childcare kit - pouch Ref 8195

P.E. Teacher first aid kit
Ref 8198



Childcare first aid kit

Ref 8170




School trip travel kit

Ref 8173

The plaster range that's safe for childminders, schools and nurseries to use on children.



There are no rules that state a responsible adult can't put a plaster on a child's minor cut. This is another example of how "Health & Safety" has been misquoted.



Ref 7184/PH

Sterile, hypo-allergenic washproof first aid adhesive dressings. Perfect for schools, could be given to all staff, (20 pack).



Rucksack first aid kit

You can order kits individually

Description	Item	Ref
Green rucksack (empty)	1	9670/GRE
Sports medium kit in rucksack	1	8090GR
Sports kit in rucksack	1	8279GR




All Steroplast's washproof hypo-allergenic plasters are perfect for a child's sensitive skin.

Some children have allergies to the adhesive used on fabric plasters (i.e. plasters that contain latex), which is why we recommend you use Steroplast's "Sterostrip" washproof plasters or our **retail range** shown above:

With so much conflicting advice on what can and can't be used on children, Steroplast have consulted with Millie's Trust, the editor of a

first aid manual and OFSTED themselves to give you a 100% guarantee that if it's got Steroplast's name on it, it is approved and safe for use on children.

The contents are based on a combination of advice from The National Childminding Association, the Health & Safety Executive, and the Pre-School Learning Alliance.



100 washproof plasters

Perfect for a child's sensitive skin supplied sterile for infection prevention

Ref 7158



4 large plasters supplied sterile

Ref 6810



First Response Kit

Grab-and-go to the emergency.

The kit enables the 'first responder' to rush to the aid of someone who's injured.

The kit opens through 180 degrees, so the contents are immediately visible, enabling an individual to administer first aid treatment faster.



New Improved Contents

Grab-and-go kit for any emergency

- The contents are clearly laid out, and secured in individual compartments when the bag is opened up fully
- Stored in water resistant bag with double zip closure - keeps contents clean and dry when used outside
- Contents held in place with elastic loops
- Adjustable nylon shoulder strap

Contents

Bag - clinical waste	1
Eyewash 500ml	1
First aid emergency leaflet	1
Foil blanket	2
Gloves - Nitrile in pairs	2
Plasters w/p	20
Resus	1
Safety Pins - pack	6
Scissors - Tufkut	1
Triangular bandage	4
Wipes - Alcohol free	10
Wound dressing - eyepad	4
Wound dressing - medium	4
Wound dressing - large	4

Description	Pack	Ref
First response kit	1	8089
First response kit - Refil	1	8089RP

40 Piece First Aid Kit

40 piece first aid kits offer unparalleled value for money.



Ref 8133

- Each compartment secured with zipper
- All products in clear compartments
- Tough and compact soft nylon case
- Comprehensive contents for general first aid

Contents	40pc
Bandage - conforming 5cm	1
Bandage - conforming 7.5cm	1
Gloves - vinyl in pairs	1
Ice pack	1
Micropore tape	1
Plasters w/p	10
Safety pins	12
Scissors - pair	1
Triangular bandage	1
Wipes - alcohol free	6
Wound dressing - pads	3
Wound dressing - adhesive	2

70 Piece First Aid Kit

Without doubt the 70 piece kits offer unparalleled quality contents.

The very same quality of products we supply into hospitals and emergency care are supplied in these comprehensive kits.

Don't take the risk of using cheaper alternatives.

- Each compartment secured with zipper
- All products are easy to see in clear compartments
- Tough and compact soft nylon roll-up case
- Velcro fastening for quick opening
- Comprehensive contents for general first aid



Ref 8128

Contents	70pc
Bandage - conforming 5cm	2
Bandage - conforming 7.5cm	2
Foil blanket	1
Gloves - vinyl in pairs	2
Ice pack	2

Contents	70pc
Micropore tape	1
Plasters w/p	30
Safety pins	12
Scissors - pair	1
Triangular bandage	2

Contents	70pc
Tweezers	1
Wipes - alcohol free	8
Wound dressing - eyepad	2
Wound dressing - pads	2
Wound dressing - adhesive	2

Fast effective treatment for burns, scalds and sunburn



Sterile Dressings

Medium Steroburn Kit
for between 1-10 people

Ref 8231

Large Steroburn Kit
for between 11-20 people

Ref 8232



Compact Steroburn Kit

Ref 8233



Sterile Dressings

Burnshield Emergency Rescue Kit

Ref 90843



Burnshield dressing in use

If your staff are split over more than one floor, you must have at least one kit per floor. First aid kits should be positioned where people see them every day, so that in the event of an emergency, people will remember where to find them.



Sterile Dressings

Mini Steroburn Kit

Ref 90842

Burns must be treated immediately!

A unique versatile foam dressing that's quick to apply - providing essential immediate treatment, relief and protection.

Remember, a burn keeps on burning until it's cooled.

Burnshield starts cooling the burn instantly. If a burn is left untreated it can lead to long term skin damage.

A one-step (remove from pack, and apply) burn treatment to reduce trauma and help prevent infection.

Burnshield hydrogel bottle

Burnshield in a bottle is ideal to treat sunburn.

Also used to rehydrate dressings that have been on the skin for an extended period of time.



125ml
Ref: BD90829



50ml Ref: BD90828

The Burnshield range of sterile dressings

Face mask	20cm x 45cm	BD90827
Finger/toes	2.5cm x 50cm	BD90825
Limb/neck	5cm x 1m	BD90826
Small area	10cm x 10cm	BD90822
Large limbs	20cm x 20cm	BD90823
Torso/back	60cm x 40cm	BD90824
Burn blotts	3.5ml x 100	BD90830
Hydrogel bottle	50ml	BD90828
Hydrogel bottle	125ml	BD90829
Suncream 1 x 10	SPF30	BD555600

Steroburn Kits

Steroburn Kit	1-10 person (M)	8231
Steroburn Kit	11-20 person (L)	8232
Steroburn Kit	Compact	8233
Steroburn Kit	Mini	90842

The Emergency Services have used Burnshield for over 20 years, giving you peace of mind that Burnshield really is the perfect product to treat skin burns effectively.

Suitable to treat 1st, 2nd or 3rd degree burns.



Suitable for warehouses, factories & kitchens

Ref 1420

Sterochef - Catering First Aid

/// The most reliable brand on the market ///



Leukoplast Coverplast X-Ray Detectable Plasters - Box of 100

The tungsten filaments contained within the plaster is its unique feature that makes it X-Ray detectable.

X-Ray detectable plasters are only used by major food manufacturers that have the very specialised equipment installed on their production lines.



100 assorted plasters in plastic case (Ref 5660)

You can rest assured that everything we supply is tried and tested.

Comes complete with a 100% satisfaction guarantee and has been designed to meet the toughest UK and European Standards.

Sterochef Blue

12 boxes of plasters/case

Description	Pack	Ref
4.0 x 4.0cm	50	5601
6.0 x 2.0cm	100	5602
7.5 x 2.5cm	100	5603
7.5 x 5.0cm	50	5605
Ass 5 sizes	100	5629
Ass 7 sizes	100	5638
Ass plastic box		5660
Fingertip	50	5912
Knuckle	50	5913
Extension extra long	50	5914
Bags of to pack into first aid kits 500/case	10	5620

Catering kits contain blue plasters

Contents	S	M	L
Bandage – conforming	1	2	2
Burn dressing	1	2	2
First aid emergency leaflet	1	1	1
Foil blanket	1	2	3
Gloves – Nitrile in pairs	6	9	12
Micropore tape	1	1	1
Plasters w/p	40	60	100
Resus device/face shield	1	1	2
Safety Pins	6	2	4
Scissors – Tuff-cut	1	1	1
Triangular bandage	2	3	4
Wall bracket	1	1	1
Wipes – saline for cleaning	20	30	40
Wound dressing – eyepad	2	3	4
Wound dressing – finger	2	3	4
Wound dressing – large	1	2	2
Wound dressing – medium	4	6	8



Sterochef selection of 7 sizes in one box - 100 plasters Ref 56385



Available in small & medium sizes

Sterochef catering kit Small Ref 8599CT-S Medium Ref 8599-CT-M

Sterochef

First aid kits, plasters and dressings for the food and catering industry.

The law requires that food processors and anyone handling food must use blue dressings.

Sterochef blue, electro-magnetically detectable dressings are coated with a medical grade, high-tack adhesive, making them the **most reliable brand on the market!**

When applied to clean dry skin, they'll stay in place, preventing them from entering the food chain, which could have very serious consequences and be very costly if a batch recall was required.

Sterochef Catering Kits, compliant with BS8599-1, supplied in a blue polypropylene 'Premier' case, or a durable blue nylon bag.

Steroplast's British Standard catering first aid kits (branded Sterochef), compliant to BS8599-1 are all you need to achieve Health & Safety compliance and peace of mind.

Description	Pack	Ref
Sterochef blue detectable plasters	1	56385
Sterochef catering kit	1	8599CT-S
	1	8599CT-M
Sterochef PE strapping	1 x 12	5682
Sterogrip blue - in sizes B, C, D, E	1 x 10	B, C, D, E
Finger stalls - Large	1 x 10	2007L
Finger stalls - Medium	1 x 10	2007M
Finger cots - flesh (latex)	1 x 10	2010L
Finger bandage - blue (see page 46)	1 x 6	3061BT



Sterogrip blue sizes Ref B, C, D, E



Sterochef strapping Ref 5682



See page 46 for blue finger bandage

Finger stalls M & L Ref 2007L

Sterogrip in blue, specifically for the catering industry. Blue Sterogrip is the essential support bandage for anyone who has sustained an injury, but continuing to work in a food and catering environment.

Sterochef Strapping is a blue plastic adhesive tape for use during food preparation, to tape around items of jewellery, preventing them from accidentally dropping into the food.

Finger Stalls & Finger Cots are easy to use and ideal for use in the food and catering industry.

The blue fingerstall must be used in food handling environments to ensure food hygiene regulation compliance.

Steroplast Premium Elastic Fabric Plaster voted #1 in the UK

Which? recommended

“Best plasters for staying power - incredibly long lasting, water resistant, sweat resistant, great value”



Ref 6038



Dressing strip 7.5cm

Ref 8002/6

Premium elastic fabric plasters and dressing strip

Our elastic fabric plaster was put through a series of tough tests by **Which?**, the 100% independent UK consumer champion to find out which plaster from major UK brands was the most reliable and cost-effective, **ours were rated #1 in the UK.**

FREE



For a complete range of retail products, order your **FREE** copy of our catalogue.

Sterostrip Washproof 12 boxes of plasters / case

Description	Pack	Ref
Spot/round	100	7165
	200	7190
4.0 x 4.0cm	50	7141
6.0 x 2.0cm	100	7142
7.5 x 2.5cm	100	7143
7.5 x 5.0cm	50	7145
Assorted 5 sizes	100	7149
Assorted 7 sizes	100	7158
Assorted plastic box	100	7160
Fingertip	50	7911
Knuckle	50	7157
Extension	50	7914
Bags of 10 to pack into first aid kits 500/case		7340
Wallet of 20 to pack into first aid kits 100/pack		7184

Premium Fabric 12 boxes of plasters / case

Description	Pack	Ref
4.0 x 4.0cm	50	6021
6.0 x 2.0cm	100	6022
7.5 x 2.5cm	100	6023
7.5 x 5.0cm	50	6025
Assorted 5 sizes	100	6059
Assorted 7 sizes	100	6038
Fingertip	50	6036
Knuckle	50	6035

*Premium Fabric Dressing Strip

Description	Pack	Ref
4.0 x 1.0m	12	6400
6.0 x 1.0m	12	6401
7.5 x 1.0m	12	6402
7.5 x 5.0m	1	6403
7.5 x 0.5m	1	8002/6

Sterochef Blue 12 boxes of plasters / case

Description	Pack	Ref
4.0 x 4.0cm	50	5601
6.0 x 2.0cm	100	5602
7.5 x 2.5cm	100	5603
7.5 x 5.0cm	50	5605
Assorted 5 sizes	100	5629
Assorted 7 sizes	100	5638
Assorted plastic box	100	5660
Fingertip	50	5912
Knuckle	50	5913
Extension	50	5914
Bags of 10 to pack into first aid kits 500/case		5620

Clear

Description	Pack	Ref
Assorted 5 sizes	100	7869

All plaster assortment boxes with tray insert, are the same size (23cm x 13cm x 6cm)



Washproof

Ref 7158



Sterocheff

Ref 5638

X-Ray detectable plasters on page 36

Blister Foam Plasters

Extra thick cushion plasters offer extra protection for treating tender areas. Ideal for corns and bunions.

- Cushioned protection - lessening pain from accidental knocks
- High-tack zinc oxide adhesive
- Amazing sticking power
- Individually wrapped - keeps plasters clean until needed
- Centralised wound pad keeps wound dry and dirt free

NEW




Ref 6050/6

There's nothing else quite like it on the market!

Our plasters stay on longer! A fact that so impressed **Which?** they voted Steroplast the No1 plaster in the UK.

1. **Lot number:** printed on the individual wrapper. This is essential for traceability, and is a **legal requirement**.
2. **Expiry date:** printed on the individual plaster. All plasters must have expiry dates printed on the reverse side of the wrapper by law!
3. **Shelf life 5 years:** from date of manufacture.
4. **Sterility:** guaranteed for 5 years
5. **CE marked:** A 100% guarantee that our products exceed the Medicines & Healthcare Regulatory Agency (M.H.R.A.) guidelines, and comply with E.U. standards.
6. **Bar coded:** all cartons are individually bar coded.
7. **Colour coded boxes:** Easy identification of box contents.
8. **Translations:** into 5 languages

Plaster sizes



- WASHPROOF**
Hypoallergenic washproof plasters
- PREMIUM**
Durable elasticated fabric plasters
- STEROFLEX**
Lightweight stretch fabric plasters
- STEROSTRIP**
Clear washproof plasters
- STEROCHEFF**
Blue detectable catering plasters

Steropore


Hypoallergenic adhesive

The most flexible adhesive wound dressing on the market.

These large sterile dressings have the benefit of a **hypoallergenic adhesive**, that has excellent adhesion, making it perfect for dressing parts of the body prone to continuous movement (e.g. joints).

Used throughout NHS hospitals following all types of surgical procedures, they are ideal because they're so flexible.

Their skin-sensitive adhesive makes them the number one choice for a child's sensitive skin, and are highly recommended for use in schools, nurseries and by childminders.



Ref 6810/1

NHS Approved



Available in 4 sizes

Ref 6806

Description	Pack	Ref
Steropore - 8.6cm x 6cm	4 x 25	6805
Steropore - 10cm x 9cm	4 x 25	6806
Steropore - 6cm x 7cm	4 x 25	6807
Steropore - 15cm x 9cm	4 x 25	6808
Steropore Plasters - 4 assorted sizes	1 x 12	6810/1

Available as a hospital or clinical pack of 4 x 25 pieces, or in a retail pack of 4 in a wallet.

Sterostrip-Sutures

Are the strongest, most reliable sutures on the market.




Ref 7858 & 7859

- The nylon strengthening fibres run the length of each suture adding to the strength of the suture, and preventing them from snapping under stress through normal body movements
- Prevent wounds from re-opening, ideal for knife cuts and similar deep wounds
- The medical grade adhesive is formulated to ensure the strip stays on longer, which in turn leads to a better healing wound or incision
- Individually wrapped and sterile



Description	Size	Pack	Ref
Sterostrip - Sutures	3mm x 75mm	10 pouches of 5	7858
Sterostrip - Sutures	6mm x 75mm	10 pouches of 3	7859



Hypoallergenic Adhesive

Sterostrip-Ultra

Ultra fine adhesive wound dressing, made from a film that conforms to the body like a second skin.

- Often referred to as a “**second-skin**” dressing because it is so thin, it’s virtually invisible
- A superb dressing chosen by surgeons following surgical procedures for improved patient comfort while allowing the individual to bath and shower

Description	Size	Pack	Ref
Sterostrip - Ultra	6cm x 7cm	1 x 25	7675/25
Sterostrip - Ultra	8cm x 10cm	1 x 25	7676/25
Sterostrip - Ultra	10cm x 12cm	1 x 25	7677/25

- Used in hospitals following surgery to cover and protect the wound site
- Sterile low adherent absorbent pad
- Allows skin to breathe
- Film coated with acrylic adhesive making it **hypoallergenic**

Steropad

Non-adhesive double-sided highly absorbent wound dressing.

Used following operations/ surgical procedures.

- The pad is double sided which takes away the risk of the pad being applied incorrectly to the wound
- Low adherent low-to-no discomfort when removed
- NHS approved
- Highly absorbent
- 4 sizes available
- Sterile



Steropad held in place with Sterotape see page 44



Steropad in boxes of 25

Ref 3007

Description	Pack	Ref
Steropad - 5cm x 5cm	4 x 25	3006
Steropad - 7.5cm x 7.5cm	4 x 25	3007
Steropad - 10cm x 10cm	4 x 25	3008
Steropad - 10cm x 20cm	25	2994BS
Dressing Pads - 5cm x 5cm	1 x 12	3006FC/1

Sterotape Premium

An unrivalled rigid tape



✓ Ultra strong

The ultra strong adhesive ensures the tape stays in place during the most rigorous sporting activity.

✓ Breathable

The specialist breathable adhesive allows the skin to breathe, making it more comfortable for the user.

✓ Ease of unwind

Unwinds with ease right down to the core, saving time & money while providing the right amount of pressure.

✓ Serrated edge

The tape tears easily and neatly across the width and equally easily down its length by hand.



DESCRIPTION	SIZE	REF
Sterotape Premium - White	2.5cm x 13.7m	4154
Sterotape Premium - Tan	3.8cm x 10m	4160/1

Find this product online or contact us to order:

Call: **0161 902 3030** | Email: sales@steroplast.co.uk

Fabric Strapping

This heavyweight tan cotton strapping is ideal for providing support whilst allowing movement. The easy to store spool and cap makes it easy to use with a smooth unwind.

SIZE	PACK	REF
2.5cm x 1.5m	1 x 12	8100FC/1



Sterotape-K Kinesiology Tape

Sterotape-K is designed to mimic the skin's elasticity. It's latex-free and low allergy.

SIZE	REF
5cm x 5m	3011 (RED)
5cm x 5m	3012 (TAN)
5cm x 5m	3010 (BLUE)



Microporous Tape

Hypo-allergenic microporous white paper tape on a simple to use spool and cap. The low allergy tape conforms to the body to securely hold dressings in place.

SIZE	PACK	REF
2.5cm x 5m	1 x 12	8105FC/1



Steroban Cohesive

The latex-free, flexible light compression bandage that sticks to itself but not the skin. It maintains consistent compression support without sticking to the patient or clothing.

SIZE	PACK	REF
4.5m	1 x 12	1630



More colours & sizes available

Sterotape Z.O.

Now with a serrated edge making it easier to tear into strips. Available in 5m and 10m lengths.

Non-stretch zinc oxide strapping

Provides rigid support even under extreme stress.

Ultra strong

Easy to apply and extremely comfortable.

Tears easily by hand along length and width

No scissors required.



Description	Pack	Ref
Sterotape-Z.O. - 1.25cm x 10m	24	4056
Sterotape-Z.O. - 2.5cm x 10m	12	4057
Sterotape-Z.O. - 4cm x 10m	8	4058
Sterotape-Z.O. - 5cm x 10m	6	4059

Leukoplast Sleek Tape

Waterproof tape

Steroplast no longer manufacture this PE tape.

So we can now offer Leukoplast waterproof tape




Description	Pack	Ref
Leukoplast Sleek tape - 2.5cm x 5m	1	TBC

Sterofix Dressing Retention Tape

Ideal for fixing dressings to the skin.

- Ideal for fixing dressings to the skin
- Non-woven polyester fabric tape
- Soft, pliable and porous
- Conforms readily to all areas of the body
- Permeable to air and water vapour

Hypoallergenic Adhesive



Description	Size	Pack	Ref
Sterofix dressing retention tape	2.5cm x 10m	4 rolls	6854
	5cm x 10m	2 rolls	6855
	10cm x 10m	1 roll	6856

Sterotape Microporous

Sterotape is made from a soft, pliable non woven material which conforms well to all contours of the body.

Spool and cap keeps tape clean and dry

Hypoallergenic Adhesive



- The adhesive layer is hypoallergenic, yet creates a secure bond when applied to the skin
- Torn easily by hand, and removed painlessly without leaving a residue
- Permeable to air and breathable

Description	Size	Pack	Ref
Microporous	1.25cm x 10m	1 x 24	6871
Microporous	2.5cm x 10m	1 x 12	6872
Microporous	5cm x 10m	1 x 6	6875
Microporous	Spool & Cap	1 x 12	8105FC-SGL

Hot & Cold Pack

Perfect for treating sports injuries

For a complete range of sports products visit www.steroplast.co.uk



Steroplast's reusable Hot and Cold pack is ideal to relieve muscular aches, pains, strains and sprains.

Orthopaedic surgeons like Steroplast's Hot and Cold pack for post-operative pain relief.

Easy to use - instructions printed on product

- One size for all - suitable for the whole family
- Reusable - a cost effective treatment for aches and pains
- Frozen pack can be used to reduce swelling and bruising on children
- Soft tactile material - nice and comfortable to use
- Washable - easy to hand wash
- Heats up easily - in your microwave or in a pan of water
- Freezes easily - put it in your freezer for approx one hour



WHEN COLD

- Reducing fevers
- For sports injuries
- Reducing swelling
- Temporary relief of headache
- Bruising
- Toothache
- Insect bites

WHEN HOT

- Temporary relief of muscle aches
- Arthritic pain
- Back pain
- Sinusitis
- Menstrual cramps
- Rheumatism

DESCRIPTION	SIZE	REF
Hot & Cold Pack	27cm x 13cm	9612/1
Hot & Cold Pack	27cm x 13cm	9612/14
Hot & Cold Pack Nylon Sleeve	27cm x 13cm	9660

Sterofreeze Ice Packs

For instant pain relief

An instant ice pack treats multiple injuries, and helps reduce swelling and bruising, and is perfect for an office, workplace or sports event as they don't need to be kept in a fridge.

No freezer or pre-cooling required

- Packs are safe to use on children
- Ideal for emergency pitch-side use
- Reduces in temperature far quicker than alternative brands
- Ideal to treat: sprains and strains immediately after the injury has happened, the most critical time
- Ice packs are ideal for all sports injuries particularly when no first aid facilities are close by
- No need to place ice packs in the freezer, simply squeezing the water sachet inside the ice pack activates them.
- Perfect for keeping in first aid kits in sports facilities and schools for unexpected injuries.



DESCRIPTION	SIZE	REF
Instant Ice Pack - Pack of 5	24cm x 12cm	9592/5
Instant Ice Pack - Case of 25	24cm x 12cm	9592
Instant Ice Pack - Pack of 10	29cm x 14cm	9630/10
Instant Ice Pack - Case of 20	29cm x 14cm	9630
Ice Pack Sleeve	-	9660

Steroply

Conforming Bandage

A high stretch, lightweight bandage for holding dressings in place.

- Used in NHS hospitals
- Soft edge for comfort
- Highly conformable, fits snugly anywhere on the body
- A nylon/viscose construction creates fray resistance
- Light and soft, permitting air circulation
- Individually wrapped
- Choice of sizes

Steroply is renowned for its quality



Description	Size	Pack	Ref
Steroply	5cm x 4m	1 x 12	3051
Steroply	7.5cm x 4m	1 x 12	3052
Steroply	10cm x 4m	1 x 12	3053
Steroply	15cm x 4m	1 x 12	3054

NHS Approved


Sterocrepe

Crepe Bandage

The NHS approved crepe bandage, ideal for giving light support to sprains and strains.

- Unrivalled quality
- Superb elasticity through "Z" and "S" twist yarns
- Can also be used for dressing retention
- Soft and supple
- Washable and therefore re-usable
- Individually wrapped

Used extensively in Health & Beauty for body wraps




Description	Size	Pack	Ref
Sterocrepe	5cm x 4.5m	1 x 12	3056
Sterocrepe	7.5cm x 4.5m	1 x 12	3057
Sterocrepe	10cm x 4.5m	1 x 12	3058
Sterocrepe	15cm x 4.5m	1 x 12	3059


Wide-open-weave (WOW) Bandage

The best selling bandages with a wide range of applications

- Loose weave construction allows wounds to breathe
- 100% cotton - Non elastic - available in two widths



Description	Size	Pack	Ref
Wow - bandage	5cm x 5m	1 x 12	3048
Wow - bandage	7.5cm x 5m	1 x 12	3049



1m tubular gauze bandage with applicator

9540/1M



Sterogauz Tubular Bandage

The ideal tubular bandage for holding dressings in place on the body.

- Ideally suited to apply over areas of the body being treated with creams and ointments, as it prevents transfer onto clothing and bedding
- Has natural stretch which provides comfort and conformability around body contours
- A range of sizes suitable for the whole body

Description	Size	Pack	Ref
Sterogauz 20m - dispensing carton	00 (1cm)	1 x 6	9545
Tubular gauze bandage	01 (1.5cm)	1 x 6	9540
	12 (2.5cm)	1 x 6	9541
	34 (3.5cm)	1	9542
	56 (5cm)	1	9543
	78 (8cm)	1	9544
Plastic applicator		25	9551/25
Plastic applicator		100	9551/100
1m tubular gauze bandage with applicator	01 (1.5cm)	1	9540/1M

Roll-on Finger Bandage

Available in blue or white

- The finger bandage can be easily applied to either fingers or toes
- The unique rolled design removes the need for applicators
- Easy to use finger bandage and cots

Description	Colour	Pack	Ref
Finger bandage	White	1 x 6	3061WT
Finger bandage	Blue	1 x 6	3061BL



1. Put roll-on finger bandage on end of finger



2. Roll the bandage down finger



3. Then twist and roll down second layer of bandage



4. Finished application



Ref 9503



Sterogrip

Elasticated Tubular Support Bandage

Elasticated tubular bandage ranging in sizes to cover ankles, knees, thighs, wrists and elbow joints.

- Designed to provide support following soft tissue injuries, dislocation and sprains. Its knitted fabric construction with interwoven elastic fibres ensures excellent support.
- The even support and pressure that Sterogrip provides make it ideal for sporting and general injuries
- Requested by the NHS specifically for its quality
- Stays in place without the need for tapes or pins
- Double layers are used for sprains, strains and over joints
- Reusable and washable
- 10 metre roll (or retail 1 metre lengths)
- Upper body sizes available on request



Retail Sterogrip 1 metre lengths

Supplied in 10m rolls or retail 1m lengths

Description	Size	Length	Ref
Sterogrip A	5cm	1 x 10m	9500
Sterogrip B	6.25cm	1 x 10m	9501
Sterogrip C	6.75cm	1 x 10m	9502
Sterogrip D	7.5cm	1 x 10m	9503
Sterogrip E	8.75cm	1 x 10m	9504
Sterogrip F	10cm	1 x 10m	9505
Sterogrip G	12cm	1 x 10m	9506
Sterogrip L	30cm	1 x 10m	9513

Sterogrip sizes


Also available in 1m lengths - minimum order 10 x 1 metre

Limb size	Child limbs	Adult limbs	Adult ankle	Adult leg	Adult thighs	Adult body
Small	A	B	C	D	E	F
Medium	B	C	D	E	F	G
Large	C	D	E	F	G	L

Sterogrip Applicator

Allows easy application over bulky dressings and double layering over painful joints.

Description	Item	Ref
Applicator (Small)	1	9537
Applicator (Medium)	1	9538
Applicator (Large)	1	9539



Prevent life threatening infections with a Sterosafe Sharps disposal kit.

Would you know how to dispose of drug related litter safely?

Used syringes, needles, blades and other sharp items can pose a serious risk to your health through infection.

- Surveys have shown that 70% of managers faced with drug related litter such as syringes, needles and blades at work would have no idea how to dispose of them safely.
- It is the manager's job to look after his/her staff by maintaining the highest standards of cleanliness and hygiene in order to be compliant with the Health & Safety at Work Act (1974).



Steroplast's Sharps disposal kit: is all you'll need to deal with drug related litter, quickly and safely. (Ref: 8205)

Provide effective handling and disposal of:

- Scalpels
- Razor blades
- Broken glass and sharp objects



Sterosafe Sharps disposal kit

- Vinyl examination gloves
- Tweezers
- Antibacterial hand wipe
- Sharp disposal - 0.08L
- Yellow disposal bag
- Steroclenz Rapid 50 spray

Ref 8205

Sterosafe Sharps disposal kit - the ultimate protection for your employees.



Wall bracket included

Ref 8410



Sharps Kits & Bins

Description	Item	Ref
Sterosafe Sharps disposal kit	1	8205
Sterosafe Sharps disposal kit	1	8410
Sharps disposal bin - 0.08L	1	8400
Sharps disposal bin - 0.3L	1	8401
Sharps disposal bin - 0.6L	1	8405
Sharps disposal bin - 1L	1	8406
Sharps disposal bin - 2L	1	8407
Sharps disposal bin - 4L	1	8409
Sharps disposal bin - 7L	1	8212

Effective handling of infectious bodily fluids with a Steroguard bio-hazard kit.



The sachet/tub of super absorbent powder is poured onto the spill **which immediately starts to turn the liquid, into a solid for easy removal.**



The yellow scoop and scraper are provided in the kit to make cleaning up the solidifying spill quick and easy. The waste can now be placed in the yellow clinical waste bag provided.



Now disinfect the area with the disinfecting spray provided. Place the used scoop, scraper and gloves in the yellow Bio-hazard bag. Dispose of yellow Bio-hazard bag according to your local authority/council guidelines.



Open the camera app on your smartphone device and hover over the QR code to find out more

Faced with highly infectious, and potentially life threatening bodily fluids to clean up, you'll need a Steroguard bio-hazard kit for your protection.

Bodily fluids including blood, vomit, urine, mucus, faeces, semen either in isolation, or on objects such as needles, syringes, broken glass and knives can all be safely disposed of using the following kits below.



Ref BIO3

Contents - Ref BIO2

- 1 **Gloves & apron:** protects against infection
- 2 **Absorbent granules:** solidifies the spill
- 3 **Scoop & scraper:** aides removal of the spill



Ref BIO1-R

- 4 **Waste bag:** for all debris
- 5 **Steroclenz 50:** for surface decontamination
- 6 **Hand wipe:** to clean hands

For safe and effective handling, containment, disinfection and disposal of infectious microorganisms contained in blood, vomit, urine, mucus, faeces and semen.

The kits meets all Health & Safety guidelines.

Steroguard bio-hazard kits

Choose 1 application for occasional use, 3 or 5 applications for frequent use.

Description	Item	Ref
Clean up kit - 5 applications	1	BIO5
Clean up kit - 3 applications	1	BIO3
Clean up pack - 1 application	1	BIO2
Clean up pack refill	1	BIO1-R
Absorbent powder - 100g tub	1	8213

First Aid Bags (Empty)



Large travel bag
24cm x 25cm

Ref 9650



Medium travel bag
20cm x 12cm oval

Ref 9651



First aid bum bag - Large
30cm x 11cm x 9cm

Ref 9653

First aid bum bag - Small
27cm x 18cm x 5cm

Ref 9659



Travel first aid belt pouch (green)
19cm x 17cm

Ref 9652



First aid belt pouch (black)
15cm x 15cm

Ref 9656



PVC first aid bag (clear green)
23cm x 18cm x 2.5cm

Ref 9664



Rucksack
27cm x 38cm x 20cm

Ref 9670/GRE



Steroplast sports bag
24cm x 29cm x 47cm

Ref DW0001



First response kit bag
35cm x 35cm x 8cm

Ref 9657

Activ Series First Aid Kit Bag

The Activ 5 Series Sports First Aid Bag is ideal for use in sports as a basic first aid kit. The bag features a carry handle and belt loop.

- Soft grip handle for comfort
- Multiple netting compartments for storage
- Belt loop on the rear for easy carriage
- Clamshell design with a double interlocking zip system



Easy to find kit contents



Belt loop on the rear



Two large zip pouches with see-through netting
Ideal to carry first aid essentials, keeping them organised in compartments.

Description	Size	Ref
Activ series 2	13.5cm x 10cm x 8cm	9672
Activ series 5	26cm x 15cm x 13cm	9675

First Aid Shoulder Bag

Ideal to carry first aid essentials, keeping them organised in compartments.

- Zipped main compartment
- Comfortable shoulder strap
- Grab handle
- Easily to identify by first aid cross
- Durable material designed to last
- Multiple compartments



Adjustable shoulder strap

Grab handle



Easy to find kit contents



Description	Size	Ref	Description	Ref
Shoulder bag (empty)	28cm x 22cm x 13cm	9654	Shoulder bag with HSE 1-10 person first aid kit contents	8018SB
			Shoulder bag with HSE 1-20 person first aid kit contents	8019SB
			Shoulder bag with BS8599 first aid kit contents	8599SB-Small
			Shoulder bag with BS8599 first aid kit contents	8599SB-Medium
			Shoulder bag with BS8599 first aid kit contents	8599SB-Large



Rucksack

A superb quality rucksack that's ideal for first aid and medical products.

Extremely versatile. Buy with contents or empty to kit out to your own specification.

Plenty of storage for extra contents: ideal for sports/schools/clubs when there is a need to have additional specific products.



Now available with first aid kit contents

- 6 inner compartments in assorted colours and sizes
- "Quick view" bags: find what you need quickly in an emergency
- Individual bags held in place with velcro for faster removal
- Document wallet on front for medical records etc
- Shoulder straps keep hands free: allows bag carrier to administer first aid
- "Sure Grip" rubber grab handle
- Grip studs on base prevents bag sitting in damp when standing upright outside
- Hi-vis reflective strips ideal for winter nights: extra safety!
- Wipe clean ultra tough material
- Dimensions approx. 33 x 39cm

'British Standard' high quality first aid kit contents, specifically designed for sports, schools, clubs and more.

Description	Ref
Rucksack with BS8599 first aid kit contents	8599RS-Small
Rucksack with BS8599 first aid kit contents	8599RS-Medium
Rucksack with BS8599 first aid kit contents	8599RS-Large
Rucksack (empty)	9668



Waterstop® Bags

Waterstop® Paramedic Emergency Backpack

The bags teflon shield has WaterStop technology to completely repel dirt and liquid.

59.5cm x 54cm x 22.5cm

- Strong carrying handles and sturdy straps
- Wipe clean PVC material
- Water-sealing zips for heavy rain
- Reflective tapes warning triangle
- Includes 4 modular bags
- Mounting for oxygen bottle



Ref PARBAG002

Waterstop® PRO Emergency Backpack (PVC or Fabric)

The bags teflon shield has WaterStop technology to completely repel dirt and liquid.

- Water-sealing zips
- Strong handles and straps
- Reflective warning triangle
- Mounting for oxygen bottle
- 2 main compartments and 2 side pockets
- Includes 5 modular bags

PVC
Ref PARBAG001

59.5cm x 54cm x 22.5cm

Fabric
Ref PARBAG001/F



Waterstop® PRO Ampoule Bag

Ampoule Holder holding up to 100 ampoules with waterproof material.

- WaterStop technology completely repels dirt and liquid
- Edge Protection

29cm x 24cm x 7.5cm



Ref AMPHOL001



Ref 9688



Steroplast Oxygen Barrel Bag

The durable wipe down oxygen barrel bag is hygienic and water-proof for all conditions

- Red oxygen barrel bag
- Internal storage pockets
- Quality fastenings
- Detachable shoulder strap

55cm x 22cm x 23cm

Steroplast Entonox Barrel Bag

High specification wipe clean PVC bag for most standard sized cylinders.

- Wipe clean PVC bag for infection control
- Strong zip opens clamshell style bag
- Interior netting ideal for storing mouthpieces and filters

55cm x 22cm x 23cm



Ref 9684

First Aid Boxes (Empty)



Holds small & medium contents - Ref 8710EB
Large kit contents - Ref 8711EB



Holds small & medium contents - Ref 8146EBC

Orange box

- Wall mountable, wall bracket supplied



Multi-purpose

Ref 8172EB

Box Dimensions:
L19 x W13.5 x H6.5cm

Premier boxes

- Durable polypropylene
- Easy access, keeps contents clean, dry and protected

Box Dimensions:

27cm x 27cm x 9cm - small & medium kit contents

34cm x 33cm x 9cm - large kit contents



Holds 1-10 & 11-20 Burns contents - Ref 8140EBC



No compartments - Clear - Ref 8152EB

Metal cabinet

- Includes two shelves
- Durable stove-enamelled construction, and lockable
- Ideal for your first aid storage



Box Dimensions:

L30cm x H46cm x W14cm

Ref 8162EBC



Free-form case with wall bracket

Box Dimensions:

26.5cm x 27.5cm x 9cm
Ref 8118EB


Wall bracket Ref 8119WB

Brackets

Wall bracket for 'Premier' boxes

Small & Medium 'Premier' boxes fit into bracket Ref 8148WB


Large 'Premier' box fits into bracket Ref 8149WB



Wall bracket
Ref 8121WB





Eye wash bracket
single bottle
Ref 8111WB



First Aid Stations

These brand new first aid stations from Steroplast come with clear printed instructions for use in an emergency.




Sterowipe

Sterile saline cleansing wipe

- Cleans skin to ensure greater adhesion of dressings
- Sterile - can be used to clean wounds
- Individual sachets (100 pack)

Description	Pack	Ref
Sterowipe - Saline cleansing wipe	1	9021



Alcohol Swabs - PDI

Sterile swabs for pre-injection, cleaning and dressing wounds

- Impregnated with a 70% Isopropyl alcohol solution
- Sterile - single use only

Description	Pack	Ref
Alcohol swabs - 30mm x 62mm	1 x 200	PDI wipes



Sterowipe

The alcohol free wipes that are ideal for cleaning the skin to ensure greater adhesion of dressings

- Supplied in a dispenser carton of 100 or in bulk packs of 2000
- Individual sachets joined in pairs
- Sterowipe is supplied in a dispenser pack for ease of use or in bulk for packing into first aid kits

Description	Pack	Ref
Sterowipe - alcohol free	1 x 100	8021
Sterowipe - alcohol free	1 x 2000	8021B



Zoff

Adhesive Remover Wipes

- Alcohol based solution for fast effective adhesive removal
- Kind to skin

Description	Pack	Ref
Zoff - adhesive remover	1 x 20	9639





Scissors & Tweezers

Ideal for a variety of medical and first aid applications.

- First aid scissors in several different styles and sizes
- Tufkut will cut quickly through tough materials and clothing in an emergency
- 7.5 inch paramedic tuff cut scissors available in 7 different colours
- Tweezers for general dressing and use by first aiders
- Allows heavily soiled dressings to be removed using tweezers rather than having direct hand contact for better hygiene
- Plastic tweezers are available sterile, and non-sterile

Description	Size	Pack	Ref
1 Nurse scissors - Blunt/Blunt	5"	1 x 10	1500
Nurse scissors - Blunt/Sharp	5"	1 x 10	1501
Nurse scissors - Sharp/Sharp	5"	1 x 10	1502
2 Medisnip - Blunt/Blunt	5"	1 x 10	1508
Medisnip - Blunt/Sharp	5"	1 x 10	1509
Medisnip - Sharp/Sharp	5"	1 x 10	1510
3 Tufkut	6"	1 x 10	1511
Tufkut	7"	1 x 10	1512
4 Megills forceps - Child	5"	1	WSTE0178
Megills forceps - Adult	10"	1	WSTE0155
Spencer Wells Artery Forceps	7"	1	WSTE0160
5 Umbilical cord scissors		1	1548



Paramedic Tufkut Scissors - Coloured

Description	Size	Item	Ref
Tufkut - Yellow scissors	7.5"	1	1512/YEL
Tufkut - Red scissors	7.5"	1	1512/RED
Tufkut - Purple scissors	7.5"	1	1512/PUR
Tufkut - Pink scissors	7.5"	1	1512/PIN

Description	Size	Item	Ref
Tufkut - Orange scissors	7.5"	1	1512/ORA
Tufkut - Blue scissors	7.5"	1	1512/BLU
Tufkut - Green scissors	7.5"	1	1512/GRE

Tweezers

Description	Size	Pack	Ref
Nickel tweezers	3"	1	1514
Martin splinter forceps	4.5"	1	1530
Plastic tweezers - Green	4.5"	10	1516
Sterile plastic tweezers - Blue	5"	10	1519/10
Sterile plastic tweezers - Blue	5"	50	1519/50

Sterile



Non-sterile



Make sure you protect yourself with PPE made to British Standards

Your best protection against viruses & bacteria

PDI Wipes

PDI Sani-Cloth Universal, the fast acting, highly effective wipes to help prevent the spread of **Coronavirus**.

Proven to kill **99.99% of germs within 30 seconds!**

Alcohol-free, PHMB free wipes for cleaning surfaces and equipment.

NHS approved

PDI wipes can be used on all types of surfaces and equipment

Extra protection against viruses & bacteria



Efficacy
Bactericidal
Yeasticidal
Virucidal (Enveloped Viruses)
Norovirus (Non-enveloped Viruses)
Tuberculocidal
Coronavirus

Description	Pack	Ref
PDI Wipes	1 x 6	PDXP00391/6

Vinyl gloves are a cost effective alternative to nitrile gloves

British Standard EN455 Parts 1, 2, 3 & 4

- Used in the food industry for hygiene reasons
- General purpose in the workplace and home



Extra protection against viruses & bacteria

- Easier on/off when dealing with customers in retail food environments

Nitrile gloves that conform to British Standards offer unrivalled protection:

British Standard EN455 Parts 1, 2, 3 & 4

- The BS EN455 series of standards highlight the tests carried out on gloves.
- Tests that cover visual defects, strength tests, checks for perforations which would render the glover ineffective at preventing bacterial transmission.



Call for more information or visit our website www.steroplast.co.uk

Nitrile powder free examination gloves

- 100 gloves per box
- Offers extra protection against certain chemicals
- Soft glove with a flexible feel
- Contains no natural rubber latex

Extra protection against viruses & bacteria

Description	Size	Pack	Ref
Nitrile blue gloves powder free	Small	1 x 100	1995
	Medium	1 x 100	1996
	Large	1 x 100	1997
	Extra Large	1 x 100	1998
Glove rack		1	1045



PDI Sani-Cloth

Recommended for the **disinfection** of non-porous hard surfaces and non-invasive medical devices

Pre-dosed disinfectant alcohol wipes providing effective infection control

Reduces the risk of infection caused by reusable cloths and liquids

Single use and ready to use



Effective against viruses & bacteria

NHS approved

Contains 70% alcohol

CALL
for
Availability

Description	Pack	Ref
70% Wipes 200 canister	1 x 6	PDXP00159/6




Aprons, coveralls & lab coat

Aprons prevents contact with an infectious agent or body fluid that may contain an infectious agent, by creating a barrier between the potential infectious material and the healthcare worker.

Description	Pack	Ref
White apron (bagged)	100	1970/100
Tyvek coverall	1	1949
Lab coat	10	1977/10

Overshoes & oversleeves

Create a barrier between the potential infectious material.



Description	Pack	Ref
Blue overshoes	1 x 200	1947
Oversleeves	1 x 100	WSTE0228

Milton sterilising fluid

- Kills 99.9% of germs including bacteria, fungi and viruses
- Sterilising solution works for over 24 hours
- Effective on Rotavirus, MRSA and swine flu virus

Extra protection against viruses & bacteria



Description	Pack	Ref
Milton 500ml	1 x 12	1430/12
Milton 1000ml	1 x 6	1433/6



Milton Sterilising Fluid is used by parents and midwives for over 70 years



Extra protection against viruses & bacteria

Hospec neutral liquid detergent

NHS
Approved

- Damp, spot and floor mopping of floors
- General cleaning and damp dusting of surfaces, hospital furniture and equipment
- Manual cleaning of dental instruments
- Bottle size is 740ml

Description	Pack	Ref
Hospec detergent	9 x 740ml	1435



First Aid Room Supplies & Essentials



First Aid Armbands (G)9677 (Y)9676



Bandage - Triangular 9562



Blanket - Red 1806



Blanket - Foil 1805



Blanket - Fleece MW1002



Torches - Pen/Pupil 1074



Ring Cutter 1290



Fob Watch WSTE0124



Whistle x2 Pack 1064

Pharmaceuticals

See page 7 for more details...



Cetrimide GSL0738



Calcium Gluconate 1959



Adhesive Label LABEL324



O.R.S. WSTE0510-13



Petroleum Jelly 9604



Security Seal 1221RE



Evac Chair WSTE0128

Disposables - Single use



Bed Pan WSTE0018/100



Bed Sheet WSTE0244



Pillowcases WSTE0081/50



Wipeable Pillow WSTE0219/10



Incontinence Pads WSTE0039



Couch Rolls 1241



Paper Towels WSTE0020



ID Bracelet 1072/100



Bic Razors 9609



Scalpels 1522



Cotton Wool Buds 2998/200



Cotton Wool Balls 3001/100



Kidney Dishes WSTE0023/100



Medicine Pots 1309/100



Medicine Spoons 1147



Tongue Depressors 1895



Vomit EM BAG 1148



Vomit Bags 1149



Vomit Bowls WSTE0037/100



Procedure Packs 1893



Sterospport **Activ** Specialist first aid for cyclists

The exceptional quality **Activ** bag ensures contents are easily identifiable, accessible and organised.

- The **Activ** range offers a comprehensive choice of first aid kits for cyclists
- Ideal kit to include in a backpack or larger bike bag

Sterospport **Activ** Adventure



Perfect for longer and more challenging cycling and outdoor adventures

Sterospport **Activ** Commuter



First aid essentials for commuter cyclists

- The excellent **Activ** bag includes everything to be prepared for longer rides
- This kit has the essentials for treating injuries experienced whilst using a bike.



Description	Reference
Sterospport Activ Essentials	SSA0001
Sterospport Activ Commuter	SSA0002
Sterospport Activ Explorer	SSA0003
Sterospport Activ Light	SSA0004
Sterospport Activ Adventure	SSA0005
Sterospport Activ Expedition	SSA0006



25 years
working with
physiotherapists
and sports
professionals

Sterosport brings a highly sought after source of knowledge, information, products and contacts to help like-minded people in the sporting community.

The first of its kind, the **"Injury Rehab Network"** started by Sterosport was set up on 3rd April 2019 as a way of connecting with clubs at all levels and to provide a comprehensive source of products and information. Sterosport wanted to get closer to key sports professionals in order that a collaboration could help us deliver even better products and services to the sports community in Manchester, the North West and nationwide.

Sterosport will help save you time and money.

When you're looking to save time and money buying tried and tested sports rehabilitation tapes and strappings, advanced life support and a reliable comprehensive range of medical essentials, or you're just after an expert's opinion and advice, Sterosport really can help!

Visit www.sterosport.co.uk, or contact us on info@sterosport.co.uk or call **0800 978 8301**

Hosting unmissable events in partnership with **BASRaT** (British Association of Sport Rehabilitators)



The **"Injury Rehab Network's"** reputation and the calibre of **"expert guest speakers"** has seen the events really take off. BASRaT's members, and the quality of the events gains audiences of over 500 attendees and growing. With unique, insightful and evidence-based presentations from world-class **"sports medicine"** professionals, the popularity continues to gain momentum. The beauty of these events is they can be attended in person or via Zoom and YouTube.

Sterosport's "Injury Rehab Network" set up to provide access to industry experts with unparalleled knowledge

SIGN UP TO OUR NEXT UNMISSABLE EVENTS...
Shows Ahead: Upcoming Guest Speakers (view all on right)

- Tom Parry**
Head of Professional Football, Manchester City FC
"Wonders of Premier League Football"
19th July 2022 | 7:00pm
- Steve Kemp**
Lead Sports Physiotherapist England Soccer Football Team
"Performance Optimisation in International Footballers"
19th July 2022 | 7:00pm
- Dr Carl Todd**
Consultant Orthopaedic Surgeon at Exeter Medical Centre
"Management of Acute Ankle Sprains"
19th July 2022 | 7:00pm

Our Previous Expert Guest Speakers Include:

- David Payne**
Head of Clinical Performance at Manchester United physio
"Team Preparation in Professional Sport"
The importance of Clinical Skills in the Modern World
No time to be equal - Pre-pandemic and Post-pandemic
- Diane Ryding**
Lead Physiotherapist for the Foundation and Youth Development Phases at Manchester United Football Club
"Building Resilient Athletes"
Hydration in an elite football academy: beyond injuries!
- Gary Lewis**
The Great Clinic - an England and Arsenal FC physio
"Rehabilitation on the road in elite sport"

Find out more...
Scan the QR code or visit www.sterosport.co.uk to register your interest for the event with Steve Kemp & other guest speakers.
Alternatively, contact Andrew Wilson.
Email: andrew@sterosport.co.uk | Tel: 0800 978 8301

Just some of the organisations that put their trust in us...

Trust is earned, and we've spent over **25 years building that trust** through our parent company Steroplast Healthcare, and more recently through its sports division Sterosport. We're proud to work in partnership with and supply products to: **professional sports teams, National Governing Bodies, community sports clubs, sports physiotherapists and a range of sports and leisure facilities.**

Put us to the test...

If you're new to **Steroplast's** sports range, and want to try it out for yourself, we'd like to encourage you with a **FREE** sample pack of tapes, a try-before-you-buy tester pack.

Just call us to take advantage of our very latest special OFFERS.

Call us on **0800 978 8301.**

To find out more and to catch up on the insights of the best of the best in sport.

Visit www.sterosport.co.uk where you can register your interest in upcoming events OR contact us on info@sterosport.co.uk

If you prefer call us on **0800 978 8301.**

"Injury Rehab Network"

was established under the Sterosport brand.



A division of Steroplast Healthcare Ltd.

▶ WATCH NEW VIDEOS EVERY WEEK

Follow our social channels for video content including:

▶ ONE MINUTE / 30 SECOND DEMOS

Quick how-to style videos demonstrating how to use our best-selling products, in 60 seconds or less.

▶ HEALTH TIPS

Learn everyday facts and tips to improve your healthcare general knowledge.

▶ NEWS

Informative videos answering your questions about current healthcare topics.

Find our videos here:



Steroplast Healthcare



@Steroplast



@Steroplast



NHS APPROVED PRODUCTS

100% SATISFACTION GUARANTEED

25 YEARS RETAIL EXPERTISE

Contact us on **0161 902 3030** for further details or to see the full range, visit **www.steroplast.co.uk**

Head Office
Steroplast Healthcare Ltd
Manchester M22 4TE England

Call: +44(0) 161 902 3030
Web: www.steroplast.co.uk
Email: info@steroplast.co.uk



11/2022