



# Alcohol free Disinfection wipes

For Surfaces and Medical Devices



Effective against NOROVIRUS in **30 Seconds** Cleans and disinfects in one

## Sani-Cloth AF Universal

are disposable disinfection wipes for non-porous hard surfaces and non-invasive medical devices. Sani-Cloth AF Universal wipes can be used to pre-clean or decontaminate medical devices prior to sterilisation or high level disinfection.

	Test Method		Contact Time
<b>Bacteria</b>	EN 16615	<i>Pseudomonas aeruginosa</i>	60 seconds
	EN 16615	<i>Staphylococcus aureus</i>	60 seconds
	EN 16615	<i>Enterococcus hirae</i>	60 seconds
	EN 13727	<i>Pseudomonas aeruginosa</i>	30 seconds
	EN 13727	<i>Staphylococcus aureus</i>	30 seconds
	EN 13727	<i>Enterococcus hirae</i>	30 seconds
<b>Enveloped Virus</b>	EN 14476	HBV, HCV, HIV, H5N1, SARS, Corona	30 seconds
	EN 14476	Vaccinia Virus - Strain Elstree	30 seconds
<b>Non-enveloped Virus</b>	EN 14476	Norovirus	30 seconds
<b>Yeasts</b>	EN 16615	<i>Candida albicans</i>	60 seconds
	EN 13624	<i>Candida albicans</i>	30 seconds
<b>TB</b>	EN 14348	<i>Mycobacterium terrae</i>	2 minutes

- Multipurpose  
Cleans and disinfects in one
- Broad spectrum of kill with low contact times
- Cost effective via product standardisation
- EN 16615**  
Supported by extracted liquid testing
- Alcohol **Free**  
PHMB **Free**
- Sani-Care**  
Training Package\*  
\*Terms & Conditions Apply



Class IIa  
Medical  
Device

200 Count Flow Wrap **Product code:** XP00390  
**Wipe size:** 267mm x 210mm

40 Count Flow Wrap **Product code:** XP00391  
**Wipe size:** 200mm x 200mm

100 Count Flow Wrap **Product code:** XP00363  
**Wipe size:** 280mm x 200mm

200 Count Canister **Product code:** XP00366  
**Wipe size:** 200mm x 200mm

225 Count Bucket **Product code:** XP00365  
**Wipe size:** 270mm x 267mm



To order call **0161 902 3030** email [marketing@steroplast.co.uk](mailto:marketing@steroplast.co.uk) or visit [www.steroplast.co.uk](http://www.steroplast.co.uk)

Setting Standards Across Seven Sectors **Emergency** Hospital Sport Veterinary Medical First Aid Retail