

# AED Plus



For more than 30 years, ZOLL® has designed products that support first responders, medical personnel and everyday rescuers in times of crisis. The ZOLL AED Plus® defibrillator, with clinically proven, guideline-driven, real-time CPR feedback and industry-leading long-life consumables, sets the benchmark for dependability, durability, performance and readiness. Bystanders will be able to act quickly and confidently in a sudden cardiac arrest emergency.



Category	Technology	ZOLL AED Plus®	Mindray BeneHeart C1A
AED	Integrated real-time CPR feedback	Yes, Real CPR Help® with: <ul style="list-style-type: none"> <li>• Compression depth bar gauge</li> <li>• Audio and visual compression depth feedback</li> <li>• Adaptive metronome</li> </ul>	No integrated, real-time CPR feedback
	Instructive prompts	<ul style="list-style-type: none"> <li>• Integrated instructive voice and text prompts</li> <li>• Visual icons guide responders through the rescue</li> </ul>	Instructive voice prompts and CPR metronome
	Pre-shock pause	13 seconds	<8 seconds
	Display (format/type)	Monochrome LCD	None
	Paediatric capabilities	Yes, intelligent paediatric capability	Yes, integrated paediatric button
	Paediatric-specific ECG analysis algorithm	Yes	No
	Paediatric-specific energy delivery	Energy levels: 50, 70, 85 J	Energy levels: 50, 70, 100 J
	Electrodes	CPR-D-padz® adult electrodes Pedi-padz® II paediatric electrodes	Universal adult/paediatric electrodes
	Electrodes: Paediatric weight/age range	<25kg or <8 years old	<25 kg or <8 years old
	Optional ECG trace	Yes	No
	Data recording and storage	1 case, 7 hours of ECG and CPR data	5 hours ECG and CPR data
	Audio recording	1 case, 20 minutes voice, ECG and CPR data	Yes, up to 1 hour
	Device data export/import	IrDA	Micro USB, USB flash memory (export), WiFi/4G (import/export)
	FDA Pre-market approval (PMA)	Yes	No

Category	Technology	ZOLL AED Plus®	Mindray BeneHeart C1A
Readiness	Self-test frequency	Weekly (default), configurable frequency for every 1, 2, 3, 4, 5, 6, and 7 days	Daily, weekly, monthly and quarterly
Managing Device Readiness	Device management software (cloud-based)	PlusTrac®	AED-Alert 2.0
	iOS/Android app for device checks	CheckAED™ mobile app	Yes
	Consumable tracking	Yes	Yes
	Responder/training tracking	Yes	No
Clinical Data Management	Case management software	<ul style="list-style-type: none"> <li>RescueNet® Event Summary</li> <li>RescueNet® CaseReview Premium</li> </ul>	No
	CPR performance summary	Yes	No
Consumables	Electrodes shelf life	<ul style="list-style-type: none"> <li>CPR-D-padz: 5 years</li> <li>Pedi-padz II: 2 years</li> </ul>	Adult and paediatric electrodes: 5 years
	Battery standby (installed)	5 years standby life	Up to 5 years standby life
	Consumables required for 10 years of use	1 set CPR-D-padz and battery	<ul style="list-style-type: none"> <li>1 set of universal electrodes</li> <li>1 battery</li> </ul>