



# Veterinary

## Essentials



# We've been the first choice of vets

for over 20 years...

The Steroplast brand may be the veterinary world's best-kept secret! We've been supplying wound care and rehabilitation products as an own label for over 20 years, so chances are you'll have used our products at some point in the recent past.



Trusted by  
vets



Manufactured to  
British Standards



Ready to use  
products



Cost-effective

## If you're a buyer, Steroplast will make your life easier

Steroplast aims to make your life as a purchaser for veterinary wound treatments and rehabilitation as easy as possible. Our products are backed up with a **100% satisfaction guarantee**, eliminating any risk associated with cheaper, inferior and unproven products.



## Vets put their trust in Steroplast for our high-quality, tried & tested:



Wound care



Hygiene & sterilisation



Clinic essentials



Tapes & bandages

Steroplast's range of adhesive dressings, tapes and bandages have also been in use throughout the NHS for more than two decades, which we think you'll agree adds further credibility to their quality and reliability.

The Steroplast brand is endorsed by healthcare professionals in all kinds of medical environments. These products are suitable to stock your veterinary practice to provide the best professional care to injured animals.

Find our full range of veterinary essentials online or contact us directly to order:

Call: **0161 902 3030** | Email: [sales@steroplast.co.uk](mailto:sales@steroplast.co.uk)

Web: [www.steroplast.co.uk](http://www.steroplast.co.uk)



## Why choose Steroplast Healthcare?

You can rest assured that everything we supply is of the highest quality, is tried and tested, comes complete with robust guarantees, gets delivered on time, and is designed to meet the toughest UK, EU, and global standards.

As a family business, we go the extra mile! Let us prove it to you.



We know there's a lot of choice in the market, but not when it comes to high-end quality, so here's our proposal...



Let us send you **two free samples** of our premium quality Elasticated Adhesive Bandage. Our EAB is a favourite with vets and in constant use across the UK & Ireland.

Based on feedback and the demand we experience, we think you'll want to try the Steroplast brand for yourself.

Please call us on **+44 (0)161 902 3030** or email [sales@steroplast.co.uk](mailto:sales@steroplast.co.uk) to claim your free sample.

# Dressings & wound care

The Steroplast ranges of adhesive dressings, tapes & bandages have been in use by the NHS for more than two decades. We think you'll agree this adds further credibility to their quality and reliability.



## ✓ High-quality

High-quality solutions for treating cuts, wounds, abrasions and lacerations.

## ✓ Ideal for emergencies

Ready to use items that veterinary professionals can use to manage major wounds in emergencies.

## ✓ Various sizes

Our wound dressings come in various sizes to reduce waste and ensure precision when treating injuries.

### Steroswab - Non-sterile

Non-sterile Steroswab is used as an additional absorbent layer for bodily fluids.

SIZE	REF
5 x 5cm	1900
7.5 x 7.5cm	1901
10 x 10cm	1902



### Steroswab - Sterile

The absorbent & sterile swab can be used when managing minor wounds - for cleaning, drying & dressing the wound.

SIZE	REF
5 x 5cm	1903
7.5 x 7.5cm	1904
10 x 10cm	1905



### Steropad wound dressings

Non-adhesive, double-sided highly absorbent wound dressings, used following operations and surgical procedures.

SIZE	REF
5 x 5cm Pack of 100	3006



More sizes available online

### Steropore wound dressings

These large sterile dressings are made with a hypo-allergenic, excellent adhesion and a highly absorbent wound pad.

SIZE	REF
10 x 9cm Pack of 25	6806/1



More sizes available online

### Wound care procedure pack

This specifically created set of disposable medical products is ideal for treating wounds and aseptic procedures in a veterinary practice.

SIZE	REF
Single pack	1893



### Sterowash eyecare kits

Our eyecare kits are perfect for flushing out any irritants or cleaning eyes or wounds.

SIZE	REF
Eye & wound	9115
Eyecare - Medium	9114
Eyecare - Large	9116



# Steroban-Supersoft

Steroban-Supersoft is a soft synthetic orthopaedic padding used to provide cushioning and extra comfort under plaster casts and as a base layer when used with pressure dressings.

## ✓ Medical Grade

Manufactured from a viscose and polyester blend, it provides a soft layer between the wound and outer dressing which helps with the healing process. The construction of the padding allows excess moisture to pass through it and evaporate, which reduces skin irritation.



PRODUCT DESCRIPTION	QTY	REF
Steroban-Supersoft 10cm x 2.7mts	12	VET003
Steroban-Supersoft 15cm x 2.7mts	12	VET004

Find this product online or contact us to order:

Call: **0161 902 3030** | Email: [sales@steroplast.co.uk](mailto:sales@steroplast.co.uk)  
Web: [www.steroplast.co.uk](http://www.steroplast.co.uk)

## ✓ Tears easily

Perfect for shaping around awkward bodily contours.

## ✓ Composition

Manufactured from a super soft viscose & polyester blend.

## ✓ Conformable

Moulds easily around any part of the body & prevents rubbing.

## ✓ Hypo-allergenic

Available in two sizes for varying coverage.

## ✓ Breathable

The lightweight padding provides improved comfort.

## ✓ Two sizes

Suitable for sensitive skin for reduced skin reactions and irritation.



# EAB Elastic Adhesive Bandage



## ✓ Premium quality elastic adhesive bandage

Steroplast's Elastic Adhesive Bandage has superb adhesive qualities, which we know from the feedback we've had from various vets is second to none.

Steroplast's EAB is available as either "normal" wound with adhesive to the inside, or "reverse wound" with adhesive to the outside.

## ✓ Central guideline

The red centre guideline on the bandage helps to keep a consistent overlap during application.

## ✓ Porous adhesive

The adhesive is porous to allow the skin to breathe and for sweat to evaporate.


## ✓ Easy unwind

Wound onto a smooth bevelled core for an easy unwind.

## ✓ Individually wrapped

All EAB is supplied individually wrapped to ensure it's clean and dry before use. Supplied in packs of 30. 2.5cm rolls in a pack of 60:

PRODUCT DESCRIPTION	QTY	REF
EAB - 2.5cm x 4.5cm	60	2092W
EAB - 5cm x 4.5cm	30	2093W
EAB - 7.5cm x 4.5cm	30	2094W
EAB - 10cm x 4.5cm	30	2095W



## ✓ Porous adhesive

Distributed evenly in stripes across the outer layer of the roll, which aids easy unwinding for a controlled, even application of pressure over painful joints.

PRODUCT DESCRIPTION	QTY	REF
Reverse wound EAB - 5cm x 4.5cm	30	2040R
Reverse wound EAB - 7.5cm x 4.5cm	30	2041R
Reverse wound EAB - 10cm x 4.5cm	30	2042R

## Reverse wound EAB

## ✓ The first choice for vets

"Reverse" wound adhesive has **adhesive on the outside** - some vets have told us that they find reverse wound makes application easier, particularly on smaller animals.

The adhesive is applied in stripes for effective ventilation and is **latex free**.

Find this product online or contact us to order:

Call: **0161 902 3030**

Email: [sales@steroplast.co.uk](mailto:sales@steroplast.co.uk)

Web: [www.steroplast.co.uk](http://www.steroplast.co.uk)

# Stero-Vet NEW Veterinary First Aid Kit

**Look after your vets when they're away from the practice.** If you are a mobile vet or often work away from your practice, you should be equipped with a personal first aid kit.

## ✓ British Standard

The contents of this kit are of the highest quality & are compliant with British Standard BS8599-1.

## ✓ Contains 500ml sterile water

This ensures an injury is adequately & hygienically cleaned before a dressing is applied. Having sterile water to hand avoids any potential for infection.

## ✓ Unique burn treatment

The burn dressing contained within the kit has a unique feature. It will cool the burn without the need for running water, relieving the pain, and allowing the vet to either keep on working or to drive back to the practice. The burn dressing can be held in place with the bandage also included in the kit.



## ✓ Designed with vets in mind

Whenever a vet is away from their practice, it makes sense to prepare for accidents that can happen in places where there may be any number of hazards, e.g. in farmyards, stables, kennels etc.

PRODUCT DESCRIPTION	REF
Stero-Vet Veterinary First Aid Kit	VET005

Find all Steroplast first aid kits online at [www.steroplast.co.uk](http://www.steroplast.co.uk), or contact us to order:

Call: **0161 902 3030** | Email: [sales@steroplast.co.uk](mailto:sales@steroplast.co.uk) | Web: [www.steroplast.co.uk](http://www.steroplast.co.uk)



# Stero-paws Animal First Aid Kit

The perfect kit for providing your family pet with pre-vet first aid treatment. It contains pet-friendly solutions, which allows pet owners to provide life-saving first aid treatments to injured animals.

## ✓ Small and lightweight

Easy to store and carry when out and about with your pet, if an emergency occurs outside of the home.

## ✓ Includes Sterowash 2-in-1

This saline solution can be used to cleanse wounds or flush out foreign objects from a pet's eyes.

PRODUCT DESCRIPTION	REF
Stero-paws Animal First Aid Kit	VET006

These kits are an ideal give-away with your **own branding**. See back page for more info on custom labels.

# Order a custom label first aid kit for your practice



## ✓ First aid kits branded with your practice logo

We can provide your practice with our high-quality first aid kits with your practice logo on them.

Get in touch today to find out how your business can take advantage of Steroplast's custom first aid kits.

## ✓ British Standard

The contents of our veterinary first aid kits are of the highest quality & are compliant with British Standard **BS8599-1**.



Contact us today to enquire about custom label kits for your own practice:

Call: **0161 902 3030** Email: [sales@steroplast.co.uk](mailto:sales@steroplast.co.uk)

**20% off**  
Steroplast®



**Special Offer:** Get **20% off** when you spend £100 or more on any Steroplast branded products.

Enter code **VET20** at checkout to apply the discount.

Plus, we'll send you a **free** printed copy of the Steroplast Healthcare Buyers Guide with your order.



Scan the QR code to shop our veterinary range, or visit [steroplast.co.uk/specialist-categories/veterinary-essentials.html](http://steroplast.co.uk/specialist-categories/veterinary-essentials.html)

## Get in touch today for more information about our veterinary supplies & own branding opportunities!

Call Steroplast today and we'll take care of everything for you. We guarantee your complete peace of mind. Starting a partnership with Steroplast is **risk free!**

Tel: **0161 902 3030** Email: [enquiries@steroplast.co.uk](mailto:enquiries@steroplast.co.uk) Web: [www.steroplast.co.uk](http://www.steroplast.co.uk)

